

Factor V Deficient Plasma immunads.

REF 5134004 Factor V Deficient Plasma, immunads.

5 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Producent / Κατασκευαστής / Производитель / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / poďtea fidiť se instrukciami / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj uputstvo pre upotrebu
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämningar / bestemmelse / Bestemmelse / προσδιορισμοί / брой тестове / stanovení / определний / ročet stanovení / Definicija
	distilled water / destilliertes Wasser / acqua distillata / água destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Serija		lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / παρτία / партида номер / šarže / lot / šarže / in vitro diagnostika
	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaksjonsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakční pufr / Reakcioni pufer		microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplatte / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrtační destička / Mikrotitrtační ploche
	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator		catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожный номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat		ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контроль / Контрольный образец / Kontrola / Kontrola		stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / späd eller upplös i / fortyndes eller opløses i / Fortyndes eller oppløses i / αρραίωση ή διάλυση σε / разтворете или разрежете с / zředí anebo rozpustí v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u		substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffer / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer		washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskoncentrat / vaskeløsningskoncentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrat promývacieho roztoku / Koncetrat solucije za ispiranje
	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico en vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рηση διαγνωστική εντόβ σωλήνα / за ин витро диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnostický prostředek in vitro / Destilšana Voda		

PRODUCT DESCRIPTION

INTENDEND USE

Factor V deficient plasma immunads. is used to determine the clotting factor V with the one stage method, also using a thromboplastin suitable for Factor V determination. The F V immune-adsorbed deficiency plasma may be used for prediluting patient specimens in determining APC resistance. The cut-off for differentiating the F V Leiden mutation has to be specifically established for the deficiency plasma – using the respective reagent and automatic equipment.

COMPOSITION

The Factor V deficient plasma immunads. is an immune-adsorbed lyophilised human deficient plasma with a clotting activity <3 %. F V Deficient Plasma immune-adsorbed is prepared from normal human plasma.

MATERIAL REQUIRED (not supplied with the kit)

-	Pipettes	-	Distilled water	-	Solutions/buffers:	
REF	5400045		Citrate-Sodium-Chloride Buffer			60 mL
-	Reagents ***					
REF	5003128		TECHNOCLLOT® PT			10 x 4 mL
REF	5003220		TECHNOCLLOT® PT Plus			5 x 4 mL
-	Calibrators and Control Plasmas ***					
REF	5220110		Coagulation Reference			5 x 1 mL
REF	5021055		Coagulation Control A			5 x 1 mL
REF	5020040		Coagulation Control N			5 x 1 mL

*** or any other package sizes, special Ceveron® alpha or Technoclot Control and Calibration reagents of Technoclone.

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- This lot of reagents prepared from human blood and each single plasma used for this lot are HBsAg, HIV 1/2 Ab and HCV Ab negative.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...+8 °C.

Stability after reconstitution:

RT* (Ceveron **)	-20 °C
2 hours	1 month

Repeated freezing is not recommended.

*= room temperature

** = in the Ceveron® alpha in the respective control area in the sample tray.

TEST PROCEDURE

PREPARATION OF PLASMA SAMPLES

Mix 9 parts of venous blood and 1 part of Sodium Citrate solution (0.11 mol/L) and centrifuge for 15 min at a RCF of at least 2500 (corresponding to DIN 58905). Store the plasma at room temperature (up to 1 day). Before carrying out the test, the patients plasma is diluted 1:10 (0.1 mL + 0.9 mL) using the Citrate-Sodium-Chloride buffer. If the presence of inhibitors is suspected, various dilutions of the plasma sample should be tested.

PERFORMANCE OF THE TEST

CEVERON

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheet supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

MANUAL

Pipetting scheme:

0.05 mL	Factor V Deficient Plasma, immune-adsorbed
+ 0.05 mL	diluted plasma samples (shake briefly and incubate for 1 min. at 37 °C)
+ 0.20 mL	thromboplastin reagent
Mix and determine the point of coagulation	

ANALYSES RESULTS

REFERENCE RANGE

70-120 % of normal

CALCULATION OF THE RESULTS

Values for 1:10 diluted plasma samples can be read off directly from the calibration curve.

When other dilutions are used, the activities in % of the calibration curve can be calculated by the following formula:

$$\frac{\% \text{ FV content}}{10} \times \text{dilution} = \% \text{ FV of the sample}$$

CALIBRATION CURVE

Reconstitute the Coagulation Reference. A predilution of 1:5 is attained by using an Citrate-Sodium-Chloride buffer in the ratio 1:5 (1 part Coagulation Reference plus 4 parts buffer). Prepare a geometric series of dilutions (1:1 to 1:64) of the predilution (1:5). The 1:1 ratio corresponds to the predilution 1:5.

Predilution	Dilutions of calibration curve						
1:5	1:1	1:2	1:4	1:8	1:16	1:32	1:64
% FV	200	100	50	25	12.5	6.3	3.1

Determine the coagulation times of the geometric series of dilutions and plot them on semilog paper (x-axis: log activity %; y-axis: coagulation time in seconds).

QUALITY CONTROL

In order to verify the accuracy of the results, Coagulation Control A and N should always be tested in the same way as the patient's plasma.

LIMITATION OF THE TEST

Incorrect sample handling can lead to partial activation of the coagulation factors and to falsely elevated single factor determinations.

Lupus anticoagulant can affect the apparent factor activity in single factor determination.

A new calibration is required for each batch of reagents where a calibration curve is necessary and for each instrument used. Also a new calibration is recommended, if software changes are introduced or following a major service of either instruments or equipment.

APPLICATIONS FOR INSTRUMENTS

Application sheets are available from Technoclone or your local distributor upon request

LITERATURE

Please contact Technoclone or see our webpage www.technoclone.com or contact your local distributor.

¹ For standardisation a reconstitution time of 30 min is recommended.

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