AA-35F1

1. Intended Use

For In Vitro Diagnostic Use

The IDS 25-Hydroxy Vitamin D RIA kit is a liquid phase radioimmunoassay intended for the quantitative determination of 25-hydroxyvitamin D (25-OH D) and other hydroxylated metabolites in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data to assist the clinician in the assessment of vitamin D sufficiency.

2. Summary and Explanation

Vitamin D is a commonly used collective term for a family of closely related seco-steroids. Upon exposure to sunlight, 7-dehydro-cholesterol, located deep in the actively growing layers of the epidermis, undergoes photolytic cleavage of the "B" ring to yield pre-vitamin D₃ which is isomerised to vitamin D₃ (cholecalciferol). Vitamin D₃ and vitamin D₂ (ergocalciferol) may also be obtained by dietary supplementation or from a limited number of foods. Vitamin D₂ is metabolised in a similar way to vitamin D₃.

Vitamin D is stored in adipose tissue and enters the circulation bound to vitamin D binding protein (VDBP) and albumin. In the liver, vitamin D is hydroxylated to give 25-hydroxyvitamin D which also circulates as a complex with VDBP. A small proportion of the 25-OH D is further hydroxylated in the kidney, under direct regulation by parathyroid hormone and ionised calcium levels, to form the biologically-active calcitropic hormone 1,25-dihydroxyvitamin D. Further hydroxylation and metabolism of vitamin D produces compounds that are water soluble and readily excreted.

Hepatic vitamin D 25-hydroxylase activity is not tightly regulated, and changes in cutaneous production of vitamin D₃, or ingestion of vitamin D (D₃ or D₂), will result in changes in circulating levels of 25-OH D ⁽¹⁾.

Serum concentration of 25-OH D is considered to be the most reliable measure of overall vitamin D status and thus can be used to determine whether a patient is vitamin D sufficient ⁽²⁾. Assessment of vitamin D status may be required to determine the cause of abnormal serum calcium concentrations in patients.

3. Method Description

The IDS 25-Hydroxy Vitamin D RIA kit is a complete assay procedure for the extraction and quantitation of 25-OH D and other hydroxylated metabolites in serum or plasma. Addition of Extraction Reagents 1 and 2 to samples and calibrators, causes precipitation of serum proteins and extraction of 25-OH D. Following centrifugation, portions of the supernatant are incubated with ¹²⁵I-labelled 25-OH D and a highly-specific sheep antibody to 25-OH D. Separation of antibody-bound tracer from free is achieved by a short incubation with Sac-Cel[®] (anti-sheep IgG cellulose) followed by centrifugation and decanting. Bound radioactivity is inversely proportional to the concentration of 25-OH D.

4. Warnings and Precautions

The IDS 25-Hydroxy Vitamin D RIA kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IDS Limited will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent. Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Radioactive Materials

- This kit contains radioactive material (Iodine-125; T_{1/2}=60.14d; γ 35.5 KeV; X 27-32 KeV). Appropriate precautions and good laboratory practices must be used in storage, handling and disposal of material. Radioactive material must be received, acquired, possessed and used only by physicians, clinical laboratories or hospitals, and only for in vitro laboratory tests. Its receipt, acquisition, possession, use and transfer are subject to local regulations.
- Store radioactive materials in the original container, in a specifically designated, properly labelled area. Access to radioactive materials must be limited to authorised personnel only.
- Do not eat, drink, smoke or apply cosmetics in areas where radioactive materials are stored or handled. Do not pipette radioactive material by mouth.
- Always wear a protective laboratory coat and disposable gloves when handling radioactive materials. Wash hands thoroughly afterwards.

- Areas where spills occur should be wiped immediately with suitable absorbent material, which should then be disposed of as radioactive waste. The contaminated area should then be washed using an alkaline detergent or radiological decontamination solution.
- Disposal of radioactive material should be in accordance with local regulations.
- Persons under 18 should not be permitted to handle radioactive material or enter radioactive areas.
- Radioactive areas must be kept clean. Use disposable or easily decontaminated laboratory ware and absorbent covers on laboratory bench surfaces to minimise contamination.

Acetonitrile

Extraction Reagent 2 **REAG** **2** contains acetonitrile.

Classification under CLP:

Acute Tox. 4

Eye Irrit. 2

Flam. Liq. 2

Hazard statements:

H225: Highly flammable liquid and vapour.

H302+312+332: Harmful if swallowed, in contact with skin or if inhaled.

H319: Causes serious eye irritation.

Precautionary statements:

P260: Do not breathe dust/fumes/gas/mist/vapours/spray.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P303+361+353: IF ON SKIN (or hair): Remove immediately all contaminated clothing. Rinse skin with water.

P301+310: IF SWALLOWED: Immediately call a POISON CENTER or doctor.

P304+340: IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

Sodium Azide

Some reagents in this kit contain sodium azide as a preservative, which may react with lead, copper or brass plumbing to form highly explosive metal azides. On disposal, flush with large volumes of water to prevent azide build up.

Sodium Hydroxide solution (<1%)

Extraction Reagent 1 **REAG** **1** contains sodium hydroxide solution (<1%).

Classification under CLP:

Eye Irrit. 2

Skin Irrit. 2

Hazard statements:

H315: Causes skin irritation.

H319: Causes serious eye irritation.

Precautionary statements:

P264: Wash hands thoroughly after handling.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332+313: If skin irritation occurs: Get medical advice/attention.

P337+313: If eye irritation persists: Get medical advice/attention.

P362: Take off contaminated clothing and wash before reuse.

Human Serum

Calibrators **CAL** and Controls **CTRL**

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV 1 and II), Hepatitis B surface antigen, antibody to Hepatitis C, and found negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled in accordance at Biosafety Level 2.

5. Shelf Life and Storage of Reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C.

Indications of possible deterioration of kit reagents

- The presence of abnormal particulate matter in any of the reagents
- A decrease in the maximum binding.
- A high non-specific binding.
- A shift in the slope of the curve from its normal position.

6. Specimen Collection and Storage

The assay should be performed using serum or plasma (EDTA or heparin) specimens. Specimens should be separated as soon as possible after collection. For long term storage, store at -20°C. Avoid repeated freeze/thaw of samples.

The specimens storage and stability information stated above are general recommendations for use in a variety of settings of laboratories. Each laboratory should follow the guidelines or requirements of local, state, and/or federal regulations or accrediting organisations to establish its own specimens handling and storage stability. For guidance on appropriate practices, please refer to the CLSI GP44-A4, Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline – Fourth Edition

7. Materials

Materials Provided

CAL 0-6 - Calibrators:	Human serum containing 25-hydroxyvitamin D and 0.09% sodium azide, 1 mL per bottle, 7 bottles per kit.
Ab - Antiserum:	Sheep anti-25-hydroxyvitamin D in phosphate buffer with 0.09% sodium azide, 52 mL per bottle, 2 bottles per kit.
Ag ¹²⁵ I - ¹²⁵ I 25-Hydroxy Vitamin D:	¹²⁵ I 25-hydroxyvitamin D in phosphate buffer with 0.09% sodium azide and proprietary stabilisers. Radioactive content <185 kBq (5 µCi) per bottle, 10.5 mL per bottle.
Sac-Cel ® :	Anti-sheep IgG coupled to cellulose suspended in buffer with 0.09% sodium azide, 11 mL per bottle.
CTRL 1 - 2 - Controls:	Human serum containing 25-hydroxyvitamin D and 0.09% sodium azide, 1 mL per bottle, 2 bottles per kit.
REAG 1 - Extraction Reagent 1:	Solution containing sodium hydroxide, <1%; 5 mL per bottle.
REAG 2 - Extraction Reagent 2:	Solution containing acetonitrile, 26 mL per bottle.
BUF - NSB Buffer:	Phosphate buffer containing 0.09% sodium azide, 8 mL per bottle.
WASHBUF 20x - Wash Concentrate:	Phosphate buffered saline containing Tween, 25 mL per bottle.

Materials Required but not Provided

- Disposable 12 x 75 mm borosilicate glass tubes.
- *Note: plastic tubes are not suitable. Do not reuse tubes.*
- Precision pipetting devices to deliver 50 µL, and 1 mL.
- Repeating pipettes to deliver 50 µL, 100 µL, 500 µL, 1 mL and 3 mL., e.g. Eppendorf Multipipette 4780, or similar.
- Vortex mixer.
- Centrifuge capable of achieving 2000 g. *A refrigerated centrifuge is preferred, but not required.*
- Decanting rack, e.g. Nalgene 5970-0013+5979-0013 or similar.
- Gamma counter capable of counting ¹²⁵I.

8. Preparation of Reagents

8.1 Allow all reagents to come to room temperature before use.

8.2 Reagents should be mixed by repeated inversion prior to use in the assay

8.3 **Wash Solution** **WASHBUF** **SOLN**

- Prepare by adding the contents of each bottle of **Wash Concentrate** **WASHBUF** **20x** to 475 mL of distilled or de-ionised water. Store at 2-8°C.

8.4 All other reagents are supplied ready for use.

9. Extraction Procedure

9.1 Prepare labelled borosilicate glass tubes, one for each Calibrator **CAL**, Control **CTRL** and unknown sample **SPE**.

9.2 Add **50 µL** of each Calibrator **CAL**, Control **CTRL** or sample to appropriately labelled tubes.

9.3 Add **50 µL** of Extraction Reagent 1 **REAG** **1** to each tube. Vortex all tubes. Incubate at 18-25°C for 20 minutes.

- 9.4 Add **500 µL** of Extraction Reagent 2 **REAG** **2** to each tube. Vortex all tubes. Centrifuge all tubes at 2000 g for 10 minutes.

10. Assay Procedure

- 10.1 Prepare labelled borosilicate glass tubes, two for each; Total Counts, NSB, Calibrator, Control and unknown sample.
- 10.2 Add **50 µL** of Calibrator, Control, or sample extract supernatants to the appropriately labelled tubes. Add 50 µL of Calibrator "A" (0) extract supernatant to NSB tubes.
- 10.3 Add **100 µL** of ¹²⁵I 25-Hydroxy Vitamin D **Ag** **¹²⁵I** to all tubes.
- 10.4 Add **1 mL** of Antiserum **Ab** to all tubes, except Total Count and NSB tubes. Add **1 mL** of NSB buffer **BUF** to NSB tubes.
- Vortex all tubes gently without foaming.
 - Incubate at 18-25°C for 90 minutes.
- 10.5 Add **100 µL** of Sac-Cel[®] to all tubes (except Total Count tubes).
- *Note: mix well to re-suspend immediately before use.*
 - Vortex all tubes gently without foaming. Incubate at 18-25°C for 20 minutes.
- 10.6 Add **3 mL** of dilute Wash Solution **WASHBUF** **SOLN** to all tubes (except Total Count tubes). Centrifuge at 2000 g for 10 minutes.
- 10.7 Decant supernatants (except Total Count tubes). Allow inverted tubes to drain on a pad of absorbent paper for 1 minute. Blot rims to remove remaining drops of liquid.
- 10.8 Count all tubes in a suitable gamma counter for at least 1 minute.

11. Calibration

25-OH D Calibrators are standardised using U.V. quantification.

12. Quality Control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Two kit controls are provided. The controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance. To ensure good precision it is important that the detectors in multi-well counters are well matched and not contaminated.

13. Calculation of Results

Calculate the percent binding (B/Bo%) of each calibrator, control and unknown sample as follows:

$$B/Bo\% = \frac{(\text{mean counts} - \text{NSB counts})}{(\text{mean counts for '0' - NSB counts})} \times 100$$

Prepare a calibration curve on semi-log graph paper by plotting B/Bo% on the ordinate against concentration of 25-hydroxy-vitamin D on the abscissa. Calculate B/Bo% for each unknown sample and read values off the curve in nmol/L (nM). Alternative data reduction techniques may be employed but users should confirm that the selected curve fit is appropriate and gives acceptable results. Smoothed spline or 4PL curve fits are recommended.

Conversion of Units:

$$\begin{array}{ccc} & \times 0.40 \Rightarrow & \\ X \text{ nmol/L} & & Y \text{ ng/mL} \\ & \Leftarrow \times 2.5 & \end{array}$$

14. Limitations of Use

- Samples suspected of containing analyte concentrations in excess of the highest calibrator should be assayed in dilution.
- Samples that contain appreciable background radioactivity should not be used. Any suspect samples should be screened for radioactivity before performing the assay. The sample should be stored until the radioactivity has decayed, or a further sample obtained.
- As in the case of any diagnostic procedure results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The following substances have been tested and found not to interfere in the 25-Hydroxy Vitamin D assay:
 - Haemoglobin: tested up to 200 mg/dL
 - Bilirubin: tested up to 500 µmol/L
 - Lipid: tested up to 3.7 mmol/L triglyceride

15. Expected Values

The following range has been determined using the IDS 25-Hydroxy Vitamin D RIA kit and is provided for guidance only. Each laboratory should determine ranges for their local population.

Normal adults 23 - 113 nmol/L (n = 24)

16. Performance Data

16.1 Accuracy

The IDS 25-Hydroxy Vitamin D RIA kit was compared against a recognised radioimmunoassay for the quantitative determination of 25-hydroxy-vitamin D and other hydroxylated metabolites. A population of 187 samples, selected to represent a wide range of 25-hydroxyvitamin D [11.8-105 nmol/L], were assayed by each method in three laboratories. Passing & Bablok regression analysis was performed on the comparative data:

IDS = $0.87(x) + 3.00$; correlation coefficient (r) = 0.95

16.2 Sensitivity

The sensitivity, defined as the concentration corresponding to the mean minus 2 standard deviations of 10 replicates of the zero calibrator, is <3 nmol/L.

16.3 Linearity

Linearity was assessed by diluting samples with buffer prior to extraction and assay.

Sample	Measured (nmol/L)	Expected (nmol/L)	% M/E
A	59.7		
A/2	27.7	28.9	92.9%
A/4	18.0	14.9	120.4%
B	43.6		
B/2	21.9	21.8	100.2%
B/4	13.0	10.9	119.1%
C	33.7		
C/2	15.2	16.8	90.3%
C/4	9.2	8.4	108.9%
		Mean	105.3%

16.4 Precision

Intra assay mean (nmol/L)	n=10 % CV	Inter assay mean (nmol/L)	n=25 % CV
26.5	5.3%	19.6	8.2%
58.4	5.0%	56.7	8.1%
151	6.1%	136	7.3%

16.5 Recovery

Recovery was assessed by adding 25-OH D₃ to samples prior to extraction and assay.

Sample conc. (nmol/L)	25-OH D ₃ added (nmol/L)	Measured (nmol/L)	Recovery (nmol/L)	Recovery %
	20	98.0	19.2	96.0%
78.8	40	119.5	40.7	101.8%
	20	51.7	19.7	98.5%
32.0	40	75.3	43.3	108.3%
	20	40.6	20.3	101.5%
20.3	40	71.4	51.1	127.8%
	20	50.2	17.8	89.0%
32.4	40	69.2	36.8	92.0%
		Mean		101.9%

16.6 Specificity

The specificity of the antiserum was assessed with the following analytes at 50% binding of the zero calibrator.

Analyte	Cross-Reactivity
25-Hydroxyvitamin D ₃	100%
25-Hydroxyvitamin D ₂	75%
24,25-Dihydroxyvitamin D ₃	>100%
Cholecalciferol D ₃	<0.01%
Ergocalciferol (D ₂)	<0.30%

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France Immunodiagnostic Systems SAS,
Germany Immunodiagnostic Systems GmbH (IDS GmbH),