DCM184-2
Ed. 04/2016

PAPP-A

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Ed. 04/2016

DCM184-2

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14

06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM184-2
Ed. 04/2016

PAPP-A

for routine analysis

Quantitative immunoenzymatic determination of PAPP-A (Pregnancy Associated Plasma Protein A) concentration in human serum.

IVD

LOT

See external label

2°C - 8°C

Σ = 96 tests

REF

DKO184

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of PAPP-A (Pregnancy Associated Plasma Protein A) concentration in human serum. PAPP-A kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

PAPP-A (Pregnancy Associated Plasma Protein A) is a 740 kDa protein produced primarily by the placenta during pregnancy; in fact this protein is usually present in the serum, but its concentration rises in maternal serum as gestation time progresses. PAPP-A concentration levels in maternal serum are indicative of threatened abortion, preterm birth, intrauterine growth restriction and others.

Many studies have shown a decrease in PAPP-A concentration in case of trisomy disorders, particularly in Down's Syndrome (Trisomy 21).

PAPP-A protein can be found in the serum in a heterotetrameric dimer with ProMBP form (most important during pregnancy) or in a dimeric form associated with coronary and cardiac conditions; Diametra PAPP-A ELISA is designed specifically for the testing of the heterotetrameric form.


When evaluation of this protein is done it is often compared to the Multiple of Medians (MoM) and represented as a percent of an in-house establish median. Without easy access to an IRP, ELISA method is an easier way to evaluate the concentration of PAPP-A in the serum and to compare the results from different laboratories.

2. PRINCIPLE

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through

the interaction of the streptavidin coated on the well and exogenously added biotinylated monoclonal anti PAPP-A antibody.

When monoclonal biotinylated antibody is mixed with a serum containing the PAPP-A antigen, an antibody-antigen complex is formed. The interaction is illustrated in the following equation:

BtAb(m) = biotinylated monoclonal antibody (excess quantity)

Ag_(PAPP-A) = native antigen (variable quantity)

Ag_(PAPP-A)-BtAb(m) = antigen-antibody complex (variable quantity)

k_a = rate constant of association

k_a = rate constant of dissociation

Simultaneously, the complex is deposited into the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

Streptavidin C.W. = Streptavidin immobilized on well
Immobilized complex = sandwich complex bound to the well

After a suitable incubation period, the antibody-antigen bound fraction is separated from unbound antigen by a washing step. Another antibody (directed at a different epitope) labelled with an enzyme is added. Another interaction occurs to form an enzyme labelled antibody-antigen-biotinylated-antibody complex on the surface of the wells.

Excess enzyme is washed off via a wash step. A suitable substrate is added to produce color measurable with the use of a microplate spectrophotometer.

The enzyme activity in the well is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can

be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATIONS

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 0.5 mL each)

CAL0	REF DCE002/18406-0
CAL1	REF DCE002/18407-0
CAL2	REF DCE002/18408-0
CAL3	REF DCE002/18409-0
CAL4	REF DCE002/18410-0
CAL5	REF DCE002/18411-0

2. Control (1 vial, 0.5 mL)

Positive serum to PAPP-A. Control concentration is indicated on the Certificate of Analysis.

REF DCE045/18403-0

3. Biotin Reagent (1 vial, 12 mL)

Biotinylated monoclonal anti PAPP-A antibody

REF DCE019/18419-0

4. Enzyme Reagent (1 vial, 12 mL)

Monoclonal anti PAPP-A antibody labeled with horseradish peroxidase (HRP)

REF DCE002/18402-0

5. Coated Microplate (1 breakable microplate)

Microplate coated with streptavidin

REF DCE002/18403-0

6. PAPP-A Diluent (1 vial, 5 mL)

Human serum based buffer

REF DCE053/18453-0

7. TMB Substrate (1 vial, 12 mL)

H₂O₂-TMB 0,26 g/L (avoid any skin contact)

REF DCE004/18404-0

8. Stop Solution (1 vial, 8 mL)

Sulphuric acid 0,15 mol/L (avoid any skin contact)

REF DCE005/18405-0

9. 50X Conc. Wash Solution (1 vial, 20 mL)

Buffered solution

REF DCE006/18406-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents between 2-8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable 2 months at 2-8°C.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.

- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the

addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The calibrators are human serum based, are obtained from pure state PAPP-A protein and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
mIU/L	0	100	500	1500	5000	10000

Once opened, the Calibrators are stable 2 months at 2-8°C.

6.2. Preparation of the Sample

The assay can be performed in human serum.

The usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain red-top venipuncture tube without additives or anti-coagulants. Allow the blood to clot, centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.020 mL of the specimen is required.

For samples with PAPP-A concentration > 10000 mIU/L, dilute the sample (for example 1:10) with PAPP-A Diluent (reagent 6); consider this dilution during the calculation of the final result.

6.3. Preparation of Wash Solution

Dilute the content of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at room temperature (22-28°C).

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each sample, one for Blank.

Reagent	Calibrator	Sample	Blank
Calibrator C ₀ -C ₅	10 µL		
Sample		10 µL	
Biotin Reagent	100 µL	100 µL	
Shake the microplate gently for 20-30 seconds to mix. Cover with a plastic wrap or microplate cover. Incubate for 60 minutes at room temperature (22-28°C). Remove the contents from each well; wash the wells 5 times with 0.350 mL of diluted wash solution. Important note: gently shake the plate for 5 seconds at each washing step to ensure proper cleaning. <u>After the last wash remove excess solution by tapping the inverted plate on an absorbent paper towel.</u>			
Enzyme Reagent	100 µL	100 µL	
Cover with a plastic wrap or microplate cover. Incubate for 30 minutes at room temperature (22-28°C). Remove the contents from each well; wash the wells 5 times with 0.350 mL of diluted wash solution. Important note: gently shake the plate for 5 seconds at each washing step to ensure proper cleaning. <u>After the last wash remove excess solution by tapping the inverted plate on an absorbent paper towel.</u>			

TMB Substrate	100 µL	100 µL	100 µL
Do not shake the plate after substrate addition. Cover with a plastic wrap or microplate cover. Incubate at room temperature (22-28°C) for 20 minutes in the dark.			
Stop Solution	50 µL	50 µL	50 µL
Shake the microplate gently for 15-20 seconds. Read the absorbance (E) at 450 nm using a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at low, normal and high range of PAPP-A for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of human PAPP-A in unknown specimens.

1. Record the absorbance obtained from the printout of the microplate reader.
2. Plot the absorbance for each duplicate serum reference versus the corresponding PAPP-A concentration in mIU/L on linear graph paper (do not average the duplicates of the serum references before plotting)
3. Draw the best fit curve through the plotted points.
4. To determine the concentration of PAPP-A for an unknown locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in mIU/L) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.586) intersects the dose response curve at 1240 mIU/L PAPP-A concentration.

Note: computer data reduction software designed for ELISA assays may also be used for the data

reduction. If such software is utilized, the validation of the software should be ascertained.

8.1. Example of Calculation

The values shown below must be considered as an example and must not be used in place of experimental data.

Calibrator /Sample	OD 450 nm	Mean OD	Value PAPP-A mIU/L
Cal 0	0.035	0.041	0
	0.046		
Cal 1	0.159	0.152	100
	0.145		
Cal 2	0.302	0.289	500
	0.276		
Cal 3	0.519	0.499	1500
	0.478		
Cal 4	1.052	1.043	5000
	1.033		
Cal 5	2.706	2.665	10000
	2.625		
Sample	0.564	0.586	1240
	0.608		

9. QUALITY CONTROL PARAMETERS

Respect the assigned ranges indicated on the Certificate of Analysis lot-specific, for both OD and Controls.

10. REFERENCE VALUES

It is recommended to compare values based on the Multiple of Median (MoM) established for the laboratory when assessing patient samples. By dividing the value of the patient sample by the MoM will give a percent value that is used frequently for evaluation.

In agreement with established reference intervals for a normal adult population, the expected ranges for Diametra PAPP-A test are:

Term of Gestation (full weeks)	PAPP-A Conc. (mIU/L)
9	650
10	1200
11	1700
12	2800
13	3900

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of

expected values based on the indigenous population where the laboratory works.

11. PERFORMANCE AND CHARACTERISTICS

11.1. Precision

The intra and inter assay precisions were determined by analyses on three different levels of pool control sera:

Intra assay				
Serum	Mean mIU/L	SD	% CV	Replicates
1	383.4	9.66	2.5	24
2	1817.2	65.69	3.6	24
3	4641.1	242.91	5.2	24

Inter assay				
Serum	Mean mIU/L	SD	% CV	Replicates
1	381.3	27.54	7.2	10
2	1900.3	97.13	5.1	10
3	4905.0	255.92	5.2	10

11.2. Sensitivity

The analytical sensitivity (detection limit) was calculated by determining the variability of the 0 mIU/L serum calibrator and using 2 standard deviation (95% certainty) statistics to calculate the minimum dose. The assay sensitivity was found to be 34 mIU/L.

11.3. Correlation

Diametra PAPP-A assay was compared with a reference immunoassay method. 50 serum samples at low, normal and high level were used.

The linear regression curve is:

$$Y = 0.5519 \cdot X + 1445.5$$

$$(R^2=0.9791)$$

11.4. Specificity

The cross-reactivity of the PAPP-A antibody to selected substances was evaluated by adding the interfering substance to a serum matrix. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of PAPP-A needed to displace the same amount of labeled analog.

Substance	Cross-reactivity
Free b-HCG	ND
AFP	ND
hPRL	ND
FSH	ND

12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Sorensen S. *Clinical Chemistry*. 2011, 57, 1023-1031.
2. Qin QP. *Clinical Chemistry*. 1997, 43, 2323-2332.
3. Fialova L., Malbohan IM. *Bratisl Lek Listy*. 2002, 103, 194-205.
4. Kagan KO. *Ultrasound Obste Gynecol*. 2008, 31, 493-502.
5. Spencer K. *Ultrasound Obste Gynecol*. 1999, 13, 231-237.
6. Baluja-Conde IB. *Rev Biomed*. 2005, 259-271.
7. Smith GCS. *The Journal of Clinical Endocrinology and Metabolism*.
8. Ong YTC. *British Journal of Obstetrics and Gynaecology*. 2000, 107, 1265-1270.
9. Bischof P. *Archives of Gynecology*. 1979, 227, 315-326

Ed. 04/2016

DCM184-2

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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