DCM181-1
Ed. 01/2015

Calcitonin

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BIBLIOGRAFIA

1. Graze K., Spiler I.J., Tashjian A.H., Melvin K.E.W., Cervi-Skinner S., Gagel R.F., Miller H.H., Wolfe H.J., Delellis R.A., Leape L., Feldman Z.T. And Reichlin S. (1978) *Natural History Of Familial Medullary Thyroid Carcinoma; Effect Of A Program For Early Diagnosis*. Engl. J. Med., 299,18;980-985.
2. Hennessy J.F., Wells S.A., Ontjes D.A. And Cooper C.W. (1974) *A Comparison Of Pentagastrin Injection And Calcium Infusion As Provocative Agents For The Detection Of Medullary Carcinoma Of The Thyroid*. J. Clin. Endocrinol. Metab., 39:487-495.
3. Rougier Ph., Calmettes C., Laplanche A., Travaglini J.P., Lefevre M., Parmentier C., Milhaud G. And Tubiana M. (1983) *The Values Of Calcitonin And Carcinoembryonic Antigen In The Treatment And Management Of Nonfamilial Medullary Thyroid Carcinoma*. Cancer, 51,5:856-862.
4. Wallach S.R., Royston I., Taetle R., Wohl H. And Deftos L. (1981) *Plasma Calcitonin As A Marker Of Disease Activity In Patients With Small Cell Carcinoma Of The Lung*. J. Clin. Endocrinol. Metab., 53,3:602-606.
5. Wells S.A., Baylin S.B., Linehan W.M., Farrell R.E., Cox E.B. And Cooper C.W. (1978) *Provocative Agents And The Diagnosis Of Medullary Carcinoma Of The Thyroid Gland*. Ann. Surg., 188,2:139-141.
6. Aurbach G.D., Marx S.J. And Spiegel A.M. (1985) *Parathyroid Hormone, Calcitonin, And The Calciferols*. In: Williams Textbook Of Endocrinology (7th Edition; Wilson J.D. And Foster D.W. Eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. Body J.J. Et Al. (1987) *ScC Antigen And Other Tumor Markers In Lung Cancer: Preliminary Results*. Excerpta Medica, 162-170.
8. Nicoli P. Et Al. (1995) *Abnormal Calcitonin Basal Levels And Pentagastrin Response In Patients With Chronic Renal Failure On Maintenance Hemodialysis*. Eur. J. Endocrinol. 132, 1, 75-81.
9. Pacini F. Et Al. (1994) *Routine Measurement Of Serum Calcitonin In Modular Disease Allows The Preoperative Diagnosis Of Unsuspected Sporadic Medullary Thyroid Carcinoma*. J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. Quesada J. M. Et Al. (1994) *Calcitriol Corrects Deficient Calcitonin Secretion In The Vit. D Deficient Elderly*. J Bone Miner Res. 1994, 9(1):53-7.

Ed. 01/2015

DCM181-1

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM181-1
Ed. 01/2015

Calcitonin

for routine analysis

Immunoenzymatic determination of Calcitonin concentration in human serum

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO181

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Calcitonin concentration in human serum.

Calcitonin kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Calcitonin is a hormone secreted by the parafollicular C cells of the thyroid gland, that acts on calcium levels by modulating the ability to "clearance" of renal calcium and reducing bone resorption. It is unclear, however, the physiological result of this modulation.

An excess of Calcitonin in the blood tissue may be indicative of the onset of malignant tumor "medullary thyroid cancer" (MTC), that mainly (80% of cases) appears sporadically, but that also has a familial component (20% of cases).

A moderate increase of Calcitonin in the blood can be found in the states of pregnancy, pernicious anemia, renal dysfunction, and during the neonatal period.

Therefore, the quantification of Calcitonin in the blood tissue can be useful mainly in the diagnosis and follow-up of MTC, but also in the study of the pathophysiology of calcium and bone metabolism.

2. PRINCIPLE OF THE METHOD

Diametra Calcitonin assay is based on the simultaneous capture of the antigen by two monoclonal antibodies (one coated on the microplate, the other linked to horseradish peroxidase) and subsequent quantification via a chromogenic substrate.

In the first phase an incubation of calibrators, controls and patient samples in the microtiter plate with the conjugate reagent is carried out; antibodies to Calcitonin form a sandwich with the antigen, blocking Calcitonin in the wells of the microplate

After incubation, the microplate is washed with a washing solution for the removal of components that have not reacted.

Finally, a chromogenic substrate solution containing TMB is added. After the incubation, the reaction is stopped by adding the Stop Solution. The color of the solution turns yellow.

The amount of color developed is directly proportional to the concentration of Calcitonin in the sample.

The concentration of Calcitonin in the sample is calculated through a calibration curve .

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, lyophilized)

CAL0	REF DCE002/18106-0
CAL1	REF DCE002/18107-0
CAL2	REF DCE002/18108-0
CAL3	REF DCE002/18109-0
CAL4	REF DCE002/18110-0
CAL5	REF DCE002/18111-0

2. Controls (1 vial, lyophilized)

Control 1	REF DCE045/18104-0
Control 2	REF DCE045/18105-0

See Control concentration on the Certificate of Analysis

3. 50X Conjugate (1 vial, 0.125 mL)

Monoclonal anti Calcitonin antibodies conjugated to horseradish peroxidase (HRP) REF DCE002/18102-0

4. Conjugate buffer (1 vial, 6 mL)

Tris buffer, BSA REF DCE002/18101-0

5. Coated Microplate (1 breakable microplate)

Microplate coated with monoclonal anti Calcitonin antibodies REF DCE002/18103-0

6. TMB Substrate (1 vial, 12 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact) REF DCE004/18104-0

7. Stop Solution (1 vial, 12 mL)

Chloridric acid 0.1N (avoid any skin contact) REF DCE005/18105-0

8. 200X Conc. Wash Solution (1 vial, 10 mL)

Tris-HCl buffer REF DCE067/18167-0

9. Serum (1 vial, lyophilized)

Human serum REF DCE053/18153-0

10. Serum buffer (1 vial, 8 mL)

Borate buffer REF DCE032/18132-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Notes

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable up to expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Thymol as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted chloridric acid solution. Chloridric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled samples.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibration (C₀...C₅)

Calibrators should be reconstituted **immediately before use with 0.5 mL of distilled or deionized water**. Shake gently until solubilization of material.

The exact concentrations of the Calibrators are stated on the labels and Certificate of Analysis for each specific lot.

Reconstituted calibrators must be stored at -20°C for a maximum of 3 months; carry out only 1 cycle of freezing and thawing (if necessary aliquot the reconstituted Calibrators).

The calibrators are calibrated against the International Standard NIBSC 89/620: 1 pg of calibrator corresponds to 0.19 µIU of International Standard.

6.2. Preparation of the Controls (1 and 2)

Proceed as the Calibrators.

Acceptability range of Controls is stated on the Certificate of Analysis.

6.3. Preparation of the Conjugate

Immediately before using, prepare the required amount of Conjugate by diluting 1:50 the Conjugate (Reagent 3) with Conjugate buffer (Reagent 4); diluted conjugate is stable 1 week at 2-8°C.

After opened, the Conjugate (reagent 3) not diluted is stable until expiry date stated on the label at 2-8°C, tightly closed.

6.4. Preparation of Wash Solution

Dilute 1:200 the "200X Conc. Wash Solution" with deionized or distilled water in the needed quantity; **prepare the diluted washing solution immediately before use and discard the diluted solution not used.**

Concentrated wash solution is stable at room temperature (22-28°C) until the expiry date stated on the label.

6.5. Preparation of the Sample

The determination of the Calcitonin with this kit can be carried out on human serum.

If the determination is not made within 24 hours of collection sample, storage the sample in aliquots at -20°C.

Avoid cycles of freezing and thawing.

Do not use hemolyzed or lipemic samples.

6.6. Preparation of the Serum (reagent 9)

Serum reagent is used to dilute any samples whose quantification is outside the calibration curve.

Reconstitute the lyophilized Serum (reagent 9) with Serum buffer (reagent 10) according to the quantity in mL indicated on the label for each specific lot.

Mix well before use.

Reconstituted Serum must be stored at -20°C for a maximum of 3 months; carry out only 1 cycle of freezing and thawing (if necessary aliquot the reconstituted Serum).

6.7. Procedure

- **Allow all reagents to reach room temperature (22-28°C).**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	100 µL		
Sample/ Control		100 µL	
Diluted Conjugate	50 µL	50 µL	
Incubate 18 (±1) hours at 2-8°C. Remove the contents from each well. Wash the wells 3 times with 400 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 30 minutes in the dark at room temperature (22÷28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Calcitonin for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (E_m) for each point of the calibration curve (C₀-C₅) and of each sample.

8.2. Calibration Curve

Plot the mean value of absorbance (E_m) of the calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.

If the sample is outside the calibration curve, dilute with reconstituted Serum (reagent 9).

9. REFERENCE VALUES

80 samples from healthy subjects were tested to establish the following normal range:

	Calcitonin (pg/mL)
Normal range	< 11 pg/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay

Within run variation was determined by replicate the measurement of two different samples in one assay. The within assay variability is $\leq 3.9\%$.

10.1.2. Inter Assay

Between run variation was determined by replicate the measurement of two different samples in different lots. The between assay variability is $\leq 6.0\%$.

10.2. Specificity

Specificity of Calcitonin kit has been checked against the following potentially interfering substances tested in concentrations of 100 ng/mL:

Substance	Interference
CGRP	No
Salmon Calcitonin	No
PDN 21	No
N-terminal procalcitonin	No

10.3. Sensitivity

The lowest detectable concentration of Calcitonin that can be distinguished from the Calibrator zero is 0.7 pg/mL at the 95% confidence limit.

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Graze K., Spiler I.J., Tashjian A.H., Melvin K.E.W., Cervi-Skinner S., Gagel R.F., Miller H.H., Wolfe H.J., Delellis R.A., Leape L., Feldman Z.T. And Reichlin S. (1978) *Natural History Of Familial*

Medullary Thyroid Carcinoma; Effect Of A Program For Early Diagnosis. Engl. J. Med., 299,18;980-985.

2. Hennessy J.F., Wells S.A., Ontjes D.A. And Cooper C.W. (1974) *A Comparison Of Pentagastrin Injection And Calcium Infusion As Provocative Agents For The Detection Of Medullary Carcinoma Of The Thyroid*. J. Clin. Endocrinol. Metab., 39:487-495.
3. Rougier Ph., Calmettes C., Laplanche A., Travagli J.P., Lefevre M., Parmentier C., Milhaud G. And Tubiana M. (1983) *The Values Of Calcitonin And Carcinoembryonic Antigen In The Treatment And Management Of Nonfamilial Medullary Thyroid Carcinoma*. Cancer, 51,5:856-862.
4. Wallach S.R., Royston I., Taetle R., Wohl H. And Deftos L. (1981) *Plasma Calcitonin As A Marker Of Disease Activity In Patients With Small Cell Carcinoma Of The Lung*. J. Clin. Endocrinol. Metab., 53,3:602-606.
5. Wells S.A., Baylin S.B., Linehan W.M., Farrell R.E., Cox E.B. And Cooper C.W. (1978) *Provocative Agents And The Diagnosis Of Medullary Carcinoma Of The Thyroid Gland*. Ann. Surg., 188,2:139-141.
6. Aurbach G.D., Marx S.J. And Spiegel A.M. (1985) *Parathyroid Hormone, Calcitonin, And The Calciferols*. In: Williams Textbook Of Endocrinology (7th Edition; Wilson J.D. And Foster D.W. Eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. Body J.J. Et Al. (1987) *ScC Antigen And Other Tumor Markers In Lung Cancer: Preliminary Results*. Excerpta Medica, 162-170.
8. Nicoli P. Et Al. (1995) *Abnormal Calcitonin Basal Levels And Pentagastrin Response In Patients With Chronic Renal Failure On Maintenance Hemodialysis*. Eur. J. Endocrinol. 132, 1, 75-81.
9. Pacini F. Et Al. (1994) *Routine Measurement Of Serum Calcitonin In Modular Disease Allows The Preoperative Diagnosis Of Unsuspected Sporadic Medullary Thyroid Carcinoma*. J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
10. Quesada J. M. Et Al. (1994) *Calcitriol Corrects Deficient Calcitonin Secretion In The Vit. D Deficient Elderly*. J Bone Miner Res. 1994, 9(1):53-7.

Ed. 01/2015

DCM181-1

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

	DE ES FR GB IT PT	In vitro Diagnostikum Producto sanitario para diagnóstico In vitro Dispositif medical de diagnostic in vitro In vitro Diagnostic Medical Device Dispositivo medico-diagnostico in vitro Dispositivos medicos de diagnostico in vitro		DE ES FR GB IT PT	Hergestellt von Elaborado por Fabriqué par Manufacturer Produttore Produzido por
	DE ES FR GB IT PT	Achtung, Begleitdokumente Precaución, consulte los documentos adjuntos Attention, veuillez consulter les documents d'accompagnement Caution, consult accompanying documents Attenzione, consultare la documentazione allegata Atenção, consultar os documentos de acompanhamento	 yyyy-mm	DE ES FR GB IT PT	Herstellungs datum Fecha de fabricacion Date de fabrication Date of manufacture Data di produzione Data de produção
 yyyy-mm-dd	DE ES FR GB IT PT	Verwendbar bis Estable hasta (usar antes de último día del mes) Utiliser avant (dernier jour du mois indiqué) Use by (last day of the month) Utilizzare prima del (ultimo giorno del mese) Utilizar (antes ultimo dia do mês)		DE ES FR GB IT PT	Biogefährdung Riesco biológico Risque biologique Biological risk Rischio biologico Risco biológico
	DE ES FR GB IT PT	Gebrauchsanweisung beachten Consultar las instrucciones Consulter le mode d'emploi Consult instructions for use Consultare le istruzioni per l'uso Consultar instruções para uso		DE ES FR GB IT PT	Chargenbezeichnung Codigo de lote Numero de lot Batch code Codice del lotto Codigo do lote
 $\Sigma = xx$	DE ES FR GB IT PT	Ausreichend für "n" Tests Contenido suficiente para "n" tests Contenu suffisant pour "n" tests Contains sufficient for "n" tests Contenuto sufficiente per "n" saggi Contém o suficiente para "n" testes		DE ES FR GB IT PT	Inhalt Contenido del estuche Contenu du coffret Contents of kit Contenuto del kit Conteúdo do kit
 Max Min	DE ES FR GB IT PT	Temperaturbereich Limitación de temperatura Limites de température de conservation Temperature limitation Limiti di temperatura Temperaturas limites de conservação		DE ES FR GB IT PT	Bestellnummer Número de catálogo Références du catalogue Catalogue number Numero di Catalogo Número do catálogo
	DE ES FR GB IT PT	Vor direkter sonneneinstrahlung schützen Mantener alejado de la luz solar Tenir à l'écart de la lumière du soleil Keep away from sunlight Tenere lontano dalla luce solare Mantenha longe da luz solar			

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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