DCM156-1
Ed. 01/2015

DIAMETRA HORMONE CONTROL PLUS

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REF DKO156

INTENDED USE

Diametra Hormone Control Plus kit provides reagents for the quality control of performance of ELISA assays in the lab, used for clinical study on patients.

1. SUMMARY AND PRINCIPLE

The use of such controls is used in the assessment of precision of method and technology used, and it is part of good laboratory practice. Three control levels are provided to guarantee a control over all the clinical range of analysis.

The product is human serum based, is lyophilized, and contains preservatives.

2. REAGENTS PROVIDED IN THE KIT

1. Controls (3 vials, lyophilized)

Control A

REF DCE045/15601-0

Control B

REF DCE045/15602-0

Control C

REF DCE045/15603-0

3. WARNINGS AND PRECAUTIONS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, Calibrators and Controls should be handled in the same manner as potentially infectious material.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- In case of microbial contamination or excessive turbidity in the reconstituted product, throw away the bottle

4. PROCEDURE

4.1. Preparation of Controls

- Allow the vials to reach room temperature (22-28°C) before using.
- Add 3 mL of deionized water to each vial. Mix well for **30 minutes**.
- Aliquot the reagent and store at -20°C.

4.2. Use of Controls

The Controls can be used on the majority of Diametra's kits (please verify the kits and the range of acceptability on the Certificate of Analysis).

For parameters such as volumes, temperature, timing, refer to the Instruction of Use of the kit.

5. STORAGE AND STABILITY

The product is stable until the expiry date stated on the label when stored at 2-8°C in the unopened original vial.

After reconstitution, all analytes within the reagent remain stable for 7 days when the reagent is stored, tightly closed, at 2-8°C, with the exception of C-peptide, Insulin and Prolactin, that remain stable only 1 day after reconstitution.

After reconstitution and subsequent freezing, all analytes are stable for 30 days if the reagent is stored tightly closed at -20°C. Once thawed, do not refreeze again but throw away what remained.

6. VALUES ASSIGNMENT

The assigned values are stated on the **Certificate of Analysis**. The values are the final result of an analysis in replicate and are **lot-specific**. The average of the obtained data must meet the range of acceptability. The average values may vary slightly depending on the population under investigation.

7. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

Ed. 01/2015

DCM156-1

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184; Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

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DiaMetra S.r.l. Headquater: Via Calabria 15,