DCM119-1
Ed. 01/2015

ENA Profile

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ENA Profile

for routine analysis


Indirect immunoenzymatic determination of IgG class autoantibodies against ENA (Extractable Nuclear Antigens) in human serum or plasma.

IVD

LOT

See external label

$\Sigma = 96$ test

REF DKO119

INTENDED USE

DiaMetra ENA Profile is an indirect solid-phase enzyme immunoassay (ELISA) for the qualitative measurement of IgG class autoantibodies against extractable nuclear antigens in human serum. The assay is intended for in vitro diagnostic use only as an aid in the diagnosis of connective tissue diseases. DiaMetra ENA Profile kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Antibodies to extractable nuclear antigens are characteristic for some rheumatic diseases and can contribute to the diagnosis and prognosis of systemic lupus erythematosus (SLE), Sjögren's syndrome, mixed-connective tissue disease (MCTD), scleroderma, poly- and dermatomyositis and CREST syndrome.

The DiaMetra ENA Profile allows for the simultaneous detection of SS-A/Ro-, SS-B/La-, SmD1-, U1-snRNP-, Scl70-, Jo1-, Histone- and Centromere B-specific autoantibodies.

SS-A / Ro Antigen: SS-A 60, SS-A 52 Sjögren's syndrome (~90%), SLE (50%)

SS-B / La Antigen: SS-B Sjögren's syndrome (85%)

SmD1 Antigen: SmD1 peptide SLE (70%)

U1-snRNP Antigen: RNP-A, RNP-C, RNP 68kD Mixed Connective Tissue Disease (100%)

Scl70 Antigen: Topoisomerase I Systemic Sclerosis (70%)

Jo1 Antigen: Histidyl-tRNA-synthetase Polymyositis / Dermatomyositis (25-30%)

Histone Antigen: Histone 2A, 2B, 3, 4 Drug-induced SLE (95 %), SLE (20-50%)

Centromere B Antigen: Centromere-associated protein B (80kD) CREST syndrome (Calcinosis, Raynaud phenomenon, Esophageal dysmotility, Sclerodactyly and Telangiectasia)

2. PRINCIPLE

The test is based on the absorptive immobilisation of extractable nuclear antigens (native, recombinant and peptidic) to microtiter strips and subsequent binding of the ENA-antibodies. The bound antibodies are detected with a peroxidase-labelled secondary antibody that is directed against human IgG. After

addition of substrate solution, a colour appears which intensity is proportional to the concentration and/or the avidity of the detected antibodies. Following the addition of stop solution, the colour switches from blue to yellow.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

- Calibrator Control (1 vial, 3 mL)
REF DCE002/11906-0
- Negative Control (1 vial, 3 mL)
Human serum based
REF DCE045/11901-0
- Conjugate (1 vial, 15 mL)
Anti human IgG conjugated with horseradish peroxidase (HRP)
REF DCE002/11902-0
- Coated Microplate (1 breakable microplate)
Coated antigens: SS-A/Ro (strip A1-12), SS-B/La (B1-12), SmD1 (C1-12), U1-snRNP (D1-12), Histone (E1-12), CENP-B (F1-12), Scl70 (G1-12) and Jo-1 (H1-12)
REF DCE002/11903-0
- Sample diluent (1 vial, 100 mL)
Phosphate buffer
REF DCE053/11953-0
- TMB Substrate (1 vial, 15 mL)
H₂O₂ 3 mM, TMB 1.2 mM (avoid any skin contact)
REF DCE004/11904-0
- Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.5M (avoid any skin contact)
REF DCE005/11905-0
- 20X Conc. Wash Solution (1 vial, 50 mL)
Tris buffer
REF DCE007/11907-0

3.2. Reagents necessary not supplied
Distilled water.

3.3. Auxiliary materials and instrumentation
Automatic dispenser.
Microplates reader (450 nm, 620-690 nm).

Note

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, the microplate is stable until the expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of standards for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, Standards and Controls should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of MIT (2-Methyl-4-isothiazolin-3-on) or Sodium Azide as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction for Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry dates printed on the labels of the box and of the vials must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user have the responsibility to make sure that the kit has been appropriately tested.

- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated should not be used in the assay. Highly lipemic or haemolysed specimens should similarly not be used
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Sample

Diametra ENA Profile assay can be performed in human serum or plasma.

Use samples freshly collected or freeze samples at -20°C. **Freeze and thaw once only.** Do not use serum samples inactivated by heat treatment at 56°C. Allow the samples to reach room temperature (30 min.).

Dilute sera 1:101 with the Sample diluent (reagent 5) (for example add 10 µL of serum to 1 mL of sample diluent)

6.2. Preparation of the Wash Solution

Dilute the content of each vial of the buffered wash solution concentrate (20X) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:20 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

6.3. Procedure

- **Allow all reagents to reach room temperature (22-28°C).**
 - Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
 - To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
1. Pipette 8 times 100 µL of diluted patient serum, the Negative Control (NC) and the Calibrator Control (CAL) in the microplate according to the following recommended scheme:

	1	2	3	...	12	antigen
A	NC	CAL	patient 1	...	patient 10	SS-A/Ro
B	NC	CAL	patient 1	...	patient 10	SS-B/La
C	NC	CAL	patient 1	...	patient 10	SmD1
D	NC	CAL	patient 1	...	patient 10	U1-snRNP
E	NC	CAL	patient 1	...	patient 10	Histone
F	NC	CAL	patient 1	...	patient 10	CENP-B
G	NC	CAL	patient 1	...	patient 10	Scl70
H	NC	CAL	patient 1	...	patient 10	Jo-1

- Seal the microplate with adhesive strip and incubate for 1h at room temperature (22-28°C)
- Discard the solution. Wash the wells 3 times with 300 µL of diluted wash solution.
- Pipette 100 µL of Conjugate and seal the microplate with adhesive strip.
- Incubate for 30 minutes at room temperature (22-28°C).
- Discard the solution. Wash the wells 3 times with 300 µL of diluted wash solution.
- Pipette 100 µL of TMB Substrate and incubate for 10 minutes at room temperature (22-28°C).
- Add 100 µL of Stop Solution.
- Read the absorbance (E) at 450 nm within the next 10 minutes after stopping. Bi-chromatic measurement with a reference wavelength at 620-690 nm is recommended.

7. QUALITY CONTROL

7.1. Validation of the test

The test results are valid provided that the following criteria are met:

- The Calibrator Control does not fall below an absorbance value of 0.3.
- The Calibrator Control does not exceed an absorbance value of 3.2
- Optionally: [CAL] > [NC]

8. RESULTS

8.1. Interpretation of Results

Interpretation of results can be made by comparing the recalculated absorbances of the Calibrator Control (CAL) and of the samples.

For each antigen the Cut-Off absorbance should be calculated according to the following equation:

OD_{Cut-Off} = OD [CAL] x Specific Factor (see Certificate of analysis)

- Absorbances > 1.1 x OD_{Cut-Off} have to be considered as positive.
- Absorbances < 0.9 x OD_{Cut-Off} have to be considered as negative.
- Absorbances 0.9 x OD_{Cut-Off} and 1.1 x OD_{Cut-Off} have to be considered as equivocal.

Calculation Example:

Row	Antigen	OD CAL	Factor	OD Cut-off	OD Sample	Result
A	SS-A/Ro	1.393	0.21	0.293	1.379	Pos
B	SS-B/La	0.562	0.55	0.309	2.442	Pos
C	SmD1	0.843	0.58	0.489	0.056	Neg
D	U1-snRNP	1.069	0.33	0.353	0.060	Neg
E	Histone	1.009	0.38	0.383	0.100	Neg
F	CENP-B	1.089	0.30	0.327	0.044	Neg
G	Scl70	1.025	0.44	0.451	0.096	Neg
H	Jo-1	0.870	0.43	0.374	0.055	Neg

9. LIMITATIONS

A positive result must be used in association with clinical evaluation and diagnostic procedures. The values obtained from this assay are intended to be an aid for diagnosis only. Elevated antinuclear antibodies may occur in individuals with no evidence of clinical disease. If the patient sample contains elevated levels of immune complexes or other immunoglobulin aggregates, false positive results by non-specific binding cannot be ruled out.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Diagnostic Sensitivity and Diagnostic Specificity

133 samples from antibody positive patients were tested to study the diagnostic sensitivity of the Diametra ENA Profile. The diagnostic specificity was determined using 100 serum samples from apparently healthy blood donors:

		Positive	Negative
Cumulative results	Pos	130	3
	Neg	3*	97
SS-A/Ro	Pos	45	0
	Neg	88	100
SS-B/La	Pos	24	0
	Neg	109	100
SmD1	Pos	18	0
	Neg	115	100
U1-snRNP	Pos	18	1
	Neg	115	99
Histones	Pos	18	2
	Neg	115	98
CENP-B	Pos	26	0
	Neg	107	100
Scl70	Pos	13	0
	Neg	120	100
Jo-1	Pos	12	0
	Neg	121	100

* Of the 3 discrepant samples testing negative in the Diametra ENA Profile but positive in ANA Screen, one sample was positive for antibodies against histones, one sample was tested positive for antibodies against SS-B/La and one sample was borderline for Scl-70 antibodies and showed a speckled, nucleolar fluorescence pattern.

Conclusion:

Relative Sensitivity	130/133	97.7%
Relative Specificity	97/100	97.0%
Relative Agreement	227/233	97.4%

10.2. Analytical Sensitivity (Detection Limit)

The analytical sensitivity (Cut-off ratio) was determined using the dilution buffer as blank. The absolute mean was obtained from 12 determinations on two different lots.

The detection limit for each antigen specificity was calculated from the absolute mean plus three standard deviations (SD). The calculated analytical sensitivities are indicated in the table below.

Antigen	Mean (OD 450nm)	SD (OD 450nm)	Analytical Sensitivity (C/O Ratio)
SS-A/Ro	0.010	0.002	0.073
SS-B/La	0.008	0.002	0.071
SmD1	0.010	0.002	0.051
U1-snRNP	0.012	0.002	0.085
Histones	0.010	0.003	0.095
CENP-B	0.010	0.003	0.084
Scl70	0.010	0.003	0.068
Jo-1	0.011	0.003	0.090

10.3. Interference

Potentially interfering substances have been evaluated by spiking three low positive samples with bilirubin at 0.2 mg/mL, hemoglobin at 400 mg/dL and lipid at 15 mg/mL.

The following table lists the tested substances in their final concentration which has been set to a level above physiologically relevant concentration.

Substance	Concentration
Bilirubin	0.2 mg/mL
Hemoglobin	4.0 mg/mL
Lipids	15 mg/mL

In general no interference could be observed by any of the tested substances with positive and negative controls up to the stated concentrations. However, in case of assaying in particular strongly hemolytic (hemoglobin > 4.0 mg/mL) and icteric (bilirubin > 0.2 mg/mL) or lipemic samples, potential interference should be taken into account leading to false results.

10.4. Precision

10.4.1. Intra-assay

For determination of the intra-assay imprecision three sera were assayed eight times each in one assay.

High positive sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	1.414	0.053	3.8
SS-B/La	0.715	0.034	4.8
SmD1	0.802	0.023	2.8
U1-snRNP	1.075	0.049	4.6
Histones	1.123	0.025	2.3
CENP-B	1.163	0.047	4.0
Scl70	1.216	0.048	3.9
Jo-1	1.103	0.021	1.9

Positive sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	1.079	0.027	2.5
SS-B/La	0.654	0.031	4.7
SmD1	0.546	0.019	3.4
U1-snRNP	0.890	0.036	4.1
Histones	0.765	0.037	4.8
CENP-B	0.793	0.022	2.8
Scl70	0.857	0.030	3.5
Jo-1	0.793	0.021	2.7

Negative sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	0.263	0.009	3.4
SS-B/La	0.094	0.005	4.9
SmD1	0.107	0.006	6.0
U1-snRNP	0.171	0.008	4.9
Histones	0.151	0.006	4.1
CENP-B	0.150	0.006	4.2
Scl70	0.136	0.004	3.2
Jo-1	0.105	0.004	3.8

The obtained data showed an intra-assay reproducibility of CV 1.9 – 6.0 %.

10.4.2. Inter-assay

For determination of the inter-assay imprecision six sera were assayed eight times each on three different runs.

High positive sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	1.269	0.06	4.9
SS-B/La	0.512	0.05	9.8
SmD1	0.703	0.03	4.0
U1-snRNP	0.930	0.04	4.6
Histones	0.897	0.02	2.2
CENP-B	0.927	0.05	5.5
Scl70	0.833	0.04	4.7
Jo-1	0.665	0.04	6.1

Positive sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	0.837	0.05	5.8
SS-B/La	0.306	0.03	8.7
SmD1	0.414	0.01	3.4
U1-snRNP	0.573	0.07	11.9
Histones	0.583	0.01	2.6
CENP-B	0.553	0.04	6.3
Scl70	0.486	0.03	6.7
Jo-1	0.387	0.03	6.7

Negative sample			
Antigen	Mean (OD 450nm)	SD	CV %
SS-A/Ro	0.248	0.01	5.8
SS-B/La	0.080	0.01	9.6
SmD1	0.099	0.01	5.2
U1-snRNP	0.162	0.02	13.4
Histones	0.148	0.01	5.4
CENP-B	0.139	0.01	9.5
Scl70	0.131	0.01	6.0
Jo-1	0.100	0.01	6.2

The obtained data confirm a good inter-assay reproducibility of CV 2.2 – 13.4%.

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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