IDK[®] DAO ELISA

***Zur in-vitro-Bestimmung von Diaminoxidase (DAO)
in Serum***

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15. LITERATUR

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). *Agents and Actions* **23**:361-65.
2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. *Scand J Clin Lab Invest* **24**:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. *Inflammation Research* **48**:169-70.

4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. *Allergy Proceedings* **15**:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. *Inflammation Research* **47**:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millenium. *Curr Probl Dermatol*, Basel, Karger, **28**:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for hiastamine induced food intolerance and supporting treatment for chronical headaches *Clin Exp Allergy* **23**: 982-85.
8. Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel *Allergologie* **9**: 426-30.
9. Jarisch Reinhart, Histamin-Intoleranz. **1. Auflage** (1999), *Thieme-Verlag*, Stuttgart

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

In-Vitro-Diagnostikum

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Hersteller

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Achtung

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Manual


IDK[®] DAO ELISA

For the in vitro determination of DAO in serum

Valid from 2018-03-22

K 8510

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE _____ 17

2. INTRODUCTION _____ 17

3. MATERIAL SUPPLIED _____ 17

4. MATERIAL REQUIRED BUT NOT SUPPLIED _____ 18

5. STORAGE AND PREPARATION OF REAGENTS _____ 19

6. SAMPLE STORAGE AND PREPARATION _____ 20

Sample stability _____ 20

Sample preparation _____ 20

7. ASSAY PROCEDURE _____ 20

Principle of the test _____ 20

Test procedure _____ 20

8. RESULTS _____ 22

9. LIMITATIONS _____ 22

10. QUALITY CONTROL _____ 23

Reference range _____ 23

Heparin treatment _____ 23

11. PERFORMANCE CHARACTERISTICS _____ 25

Sample dilution _____ 25

Analytical Sensitivity _____ 25

Precision and reproducibility _____ 26

12. PRECAUTIONS _____ 26

13. TECHNICAL HINTS _____ 27

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE _____ 27

15. REFERENCES _____ 27

1. INTENDED USE

This Immundiagnostik AG assay is an enzyme immunoassay intended for the quantitative determination of diamine oxidase (DAO) in serum. For *in vitro* diagnostic use only.

2. INTRODUCTION

Diamine oxidase (DAO) is a body's own enzyme that metabolises histamine. Although DAO is found practically in the whole body, the most important site of its action is the intestine. The enzymatic activity of DAO determines the histamine degradation speed. In the case of DAO deficiency or inhibition, incorporated or endogenous histamine cannot be degraded quickly enough, and the symptoms of histamine intolerance are presented. Millions of people suffer from gastrointestinal problems, migraine, irritations of nasal mucosa and other allergy-like symptoms after consumption of certain nutrients. Too much histamine in the body can be the reason for this wide range of symptoms.

The determination of DAO serum concentration (K 8510) combined with the determination of DAO activity (K 8220 DAO REA) is a suitable marker for the differential diagnosis of histamine intolerance and associated symptoms.

Our IDK® DAO ELISA kit is intended for determination of the diamine oxidase (DAO) concentration in serum.

Indications

- Frequent headaches or migraine
- Snuffles after consumption of histamine-containing nutrients
- Tissue oedema
- Eyelid turgor
- Skin redness
- Limb aches
- Gastrointestinal discomfort
- Monitoring of a histamine free diet

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 8510	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 8510	WASHBUF	Wash buffer concentrate, 5 x	4 x 100 ml

Cat. No.	Label	Kit components	Quantity
K 8510	CAL	Calibrator, lyophilised (see specification for concentration)	4 x 1 vials
K 8510	CTRL1	Control, lyophilised (see specification for range)	4 x 1 vial
K 8510	CTRL2	Control, lyophilised (see specification for range)	4 x 1 vial
K 8510	AB	Detection antibody concentrate , lyophilised (biotinylated)	1 x 200 µl
K 8510	CONJ	Conjugate concentrate, peroxidase-labelled	1 x 200 µl
K 8510	ABBUF	Dilution buffer for AB und CONJ, ready-to-use	1 x 50 ml
K 8510	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 50 ml
K 8510	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 8510	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Standard laboratory reaction vessels (1.5 ml)
- Standard laboratory reaction vessel (15 ml)
- Foil to cover the microtiter plate
- Microtiter plate thermoshaker at 37 °C (for example model Shake ID2 available at Immundiagnostik AG)
- Vortex mixer
- Ultra pure water*
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:5** before use (200 ml WASHBUF + 800 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:5 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

Please note:

This WASHBUF is intended only for use in the IDK® DAO ELISA. Crystals in the WASHBUF must be completely dissolved before dilution.

- The lyophilised **CAL** (calibrator) and **CTRL** (controls) are stable at 2–8 °C until the expiry date stated on the label. Before use, the CAL and CTRL have to be reconstituted with **500 µl sample dilution buffer (SAMPLEBUF)** and mixed by gentle inversion to ensure complete reconstitution. Allow the vial content to dissolve for 10 minutes and then mix thoroughly. **Calibrator and controls** (reconstituted CAL and CTRL) **are not stable and cannot be stored**.
- Use **100 µl SAMPLEBUF** (sample dilution buffer) as **BLANK**.
- **Preparation of the conjugate and the detection antibody:** Before use, the **conjugate concentrate (CONJ)** and the **detection antibody concentrate (AB)** have to be diluted **1:101** in **dilution buffer** (100 µl CONJ + 10 ml ABBUF), (100 µl AB + 10 ml ABBUF). The CONJ and the AB are stable at **2–8 °C** until the expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) and **detection antibody** (1:101 diluted AB) **are not stable and cannot be stored**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

6. SAMPLE STORAGE AND PREPARATION

Sample stability

Storage at room temperature (15–30 °C):	4 days
Storage at 2–8 °C:	9 days
Storage at -20 °C:	6 months

Sample preparation

Serum samples must be diluted **1:5** before performing the assay,

e.g. **50 µl** sample + **200 µl** sample dilution buffer (SAMPLEBUF), mix well.

100 µl of the dilution are used in the test.

7. ASSAY PROCEDURE

Principle of the test

This ELISA is designed for the quantitative determination of DAO in serum. The assay utilises the “sandwich” technique with two polyclonal antibodies against recombinant DAO.

Calibrator, controls and diluted samples which are assayed for DAO are added into the wells of a micro plate coated with polyclonal rabbit anti- DAO antibody. During the first incubation step, DAO is bound by the immobilised primary antibody. Then a biotinylated polyclonal anti-DAO antibody, is added into each microtiter well. In the next step, the streptavidin-POD-conjugate is added and a “sandwich” of

1st antibody – DAO - biotinylated antibody – streptavidin-POD-conjugate is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of DAO. The concentration of DAO can be quantified by referring the optical density of the calibrator to a lot-dependent master calibration curve.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of calibrator/controls/blank/samples on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips in the carefully closed aluminium packaging at 2–8° C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Before use , wash the wells 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Add each 100 µl calibrator/controls/blank/diluted samples into the respective wells.
3.	Cover the strips and incubate for 2 hours at 37 °C on a horizontal shaker* .
4.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5.	Add 100 µl detection antibody (diluted AB) into each well, mix gently.
6.	Cover the strips and incubate for 1 hour at 37 °C on a horizontal shaker* .
7.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8.	Add 100 µl conjugate (diluted CONJ) into each well.
9.	Cover the strips and incubate for 1 hour at 37 °C on a horizontal shaker* .
10.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
11.	Add 100 µl substrate (SUB) into each well.
12.	Incubate for 10–20 min** at room temperature (15–30 °C) in the dark .
13.	Add 100 µl stop solution (STOP) into each well and mix well.
14.	Determine absorption immediately with an ELISA reader at 450 nm . If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

* We recommend shaking the strips at 700 rpm with an orbit of 2 mm.

** The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

8. RESULTS

For result evaluation, please use a four parametric logit-log model based on the standard curve of the respective kit lot and the calibrator value (CAL). All essential information on the standard curve is provided on the QC data sheet of the respective product lot.

The calibration curve can be expressed either by the concentration of each standard with its corresponding optical density or by the four parameters A,B,C and D. In both cases the optical density of the calibrator (CAL) is essential. Depending on your evaluation software program, either the one or the other kind of data described above should be entered.

Caution: Please make sure that all parameters and values are transferred accurately into your software as minor deviations can cause severe errors during evaluation.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Samples

The obtained results have to be multiplied by the **dilution factor of 5** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the calibration curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

LoB see chapter "Performance Characteristics".

10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range


< 3 U/ml:	high incidence for HIT (Histamine intolerance)
3 - 10 U/ml:	HIT probable
> 10 U/ml:	low incidence for HIT

Conversion factor: 1 U/ml = 1.25 ng/ml

We recommend each laboratory to establish its own reference range.

Comparison “Pregnant women” and “Healthy controls”

For the clinical evaluation of this assay we have analysed samples from pregnant women and apparently healthy controls. The IDK® DAO ELISA detects, as required and expected, higher values in pregnant women than in healthy controls.

Heparin treatment

Furthermore, DAO levels in healthy study participants increased sharply within 30 minutes of heparin administration. It has been documented in scientific literature that DAO levels rise after heparin administration.

Treatment outcome before and after heparin administration

IDK® DAO ELISA [U/ml]

	BEFORE administration	30 min AFTER	60 min AFTER
Patient 1	76	219	-
Patient 2	55	152	-
Patient 3	2.5	277	622
Patient 4	18.9	621	555

Since the IDK® DAO ELISA determines DAO concentration while the conventional histamine intolerance tests with putrescine or histamine as substrates determine DAO activity, the correlation coefficient must not necessarily be $r > 0.8$. This can be explained by the fact that the activity does not depend on the number of molecules alone, but also on cofactors such as vitamin C, vitamin B6, copper or manganese ions in vitro and in vivo. For the diagnosis of histamine intolerance via DAO activity test we therefore recommend to determine the above mentioned cofactors as well. The problem may not be a low DAO level, but a cofactor deficiency.

The symptoms of histamine intolerance can be caused by low DAO activity because the above-mentioned cofactors are not sufficiently available. By quantitating the cofactors it can be determined which one needs to be supplemented.

Medication effects

In addition, histamine intolerance symptoms may be due to low DAO activity caused by medication such as:

Muscle relaxants	Pancuronium, alcuronium, D-tubocurarine
Narcotics	Thiopental
Analgetics	Morphine, pethidine, nonsteroidal anti-inflammatory drugs, acetylsalicylic acid, metamizole
Local anesthetics	Prilocaine
Antihypotonics	Dobutamine
Antihypertensive drugs	Verapamil, alprenolol, dihydralazine
Antiarrhythmics	Propafenone
Diuretics	Amiloride
Drugs influencing gut motility	Metoclopramide

Antibiotics	Cefuroxime, cefotiam, isoniazid, pentamidin, clavulanic acid, choroquine
Mucolytics	Acetylcysteine, ambroxol
Broncholytics	Aminophylline
H2-receptor antagonists	Cimetidine
Cytostatics	Cyclophosphamide
Antidepressants	Amitriptyline

If you are taking such medication, you may want to discuss with your physician alternative medication in order to relieve your symptoms.

11. PERFORMANCE CHARACTERISTICS

Sample dilution

Two patient samples were diluted and assayed. The results are shown below:

Linearity n= 2

Sample	Dilution	Expected [U/ml]	Measured [U/ml]
A	1:5	15.65	15.65
	1:10	7.83	7.68
	1:20	3.92	3.86
	1:40	1.96	1.86
B	1:5	5.58	5.58
	1:10	2.79	2.89
	1:20	1.39	1.45
	1:40	0.7	0.7

Analytical Sensitivity

Limit of blank, LoB	0.079 U/ml
Limit of detection, LoD	0.120 U/ml
Limit of quantitation, LoQ	0.120 U/ml

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20 % CV.

Precision and reproducibility

Intra-Assay (n = 16)

The precision (intra-assay variation) was calculated from 16 replicate determinations on each of two samples.

Sample	DAO [U/ml]	CV [%]
1	5.0	1.42
2	5.9	1.72

Inter-Assay (n = 8)

The total precision (inter-assay variation) was calculated from data on 2 samples obtained in 8 different assays by different technicians on three different days.

Probe	DAO [U/ml]	CV [%]
1	23.27	7.9
2	13.36	10.7

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDK® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Sattler J et al. (1988) Food induced histaminosis as an epidemiological problem: plasma histamine elevation and haemodynamic alterations after oral histamine administration and blockade of diamine oxidase (DAO). *Agents and Actions* **23**:361-65.

2. Tufvesson G et al. (1969) Determination of DAO-activity in normal human blood serum. *Scand J Clin Lab Invest* **24**:163-68.
3. Wantke F et al. (1999) The red wine maximization test: drinking histamine rich wine induces a transient increase of plasma diamine oxidase activity in healthy volunteers. *Inflammation Research* **48**:169-70.
4. Wantke F et al. (1994) The red wine provocation test: intolerance to histamine as a model für food intolerance. *Allergy Proceedings* **15**:27-32.
5. Wantke F et al. (1998) Daily variations of serum diamine oxidase and the influence of H1 and H2 blockers: a critical approach to routine diamine oxidase assessment. *Inflammation Research* **47**:396-400.
6. Jarisch R et al. (1999) Role of food allergy and food intolerance in recurrent urticaria. In: Wüthrich B (Hrsg): The Atopy Syndrome in the Third Millenium. *Curr Probl Dermatol*, Basel, Karger, **28**:64-73.
7. Wantke F et al. (1993) Histamine free diet: treatment of choice for hiastamine induced food intolerance and supporting treatment for chronical headaches *Clin Exp Allergy* **23**: 982-85.
8. Götz M et al. (1996) Histamin-Intoleranz und Diaminoxidasemangel *Allergologie* **9**: 426-30.
9. Jarisch Reinhart, Histamin-Intoleranz. **1. Auflage** (1999), *Thieme-Verlag*, Stuttgart

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by

Attention

Consult instructions for use

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49 (0) 62 51/70 19 00

Fax: +49 (0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com