CRP ELISA Kit

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15. LITERATUR

1. Koenig, W. et al., 2004. C-reactive protein modulates risk prediction based on the Framingham Score: implications for future risk assessment: results from a large cohort study in southern Germany. *Circulation*, **109**(11), pp.1349–53.
2. Pearson, T.A., 2003. American Heart Association Guide for Improving Cardiovascular Health at the Community Level: A Statement for Public Health Practitioners, Healthcare Providers, and Health Policy Makers From the American Heart Association Expert Panel on Population and Pre. *Circulation*, **107**(4), pp.645–651.
3. Ridker, P.M. et al., 2000. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. *The New England journal of medicine*, **342**(12), pp.836–43.

Verwendete Symbole:

Temperaturbegrenzung

In-Vitro-Diagnostikum

Hersteller

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Bestellnummer

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
CRP ELISA Kit

***For the in vitro determination of C-reactive protein from
dried blood samples***

Valid from 2015-10-01

K9710SDBS

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE 15

2. INTRODUCTION 15

3. MATERIAL SUPPLIED 15

4. MATERIAL REQUIRED BUT NOT SUPPLIED 16

5. STORAGE AND PREPARATION OF REAGENTS 16

6. PREPARATION OF SAMPLES 17

7. ASSAY PROCEDURE 17

Principle of the test 17

Test procedure 17

8. RESULTS 19

9. LIMITATIONS 19

10. QUALITY CONTROL 19

Reference range 20

11. PERFORMANCE CHARACTERISTICS 20

Specificity 20

Analytical Sensitivity 20

Spiking Recovery 21

Precision and reproducibility 21

Dilution recovery 22

12. PRECAUTIONS 22

13. TECHNICAL HINTS 23

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE 23

15. REFERENCES 24

1. INTENDED USE

The Immundiagnostik Assay is a sandwich enzyme immuno assay intended for the quantitative determination of C-reactive Protein in dried blood samples. It is for in vitro diagnostic use only.

2. INTRODUCTION

C-reactive protein (CRP) is mainly formed in hepatocytes. The synthesis rate of CRP is influenced by the cytokines involved in the inflammatory processes. The biological half-life time is estimated to be 13-16 hours. The serum CRP concentration reflects very sensitive acute fever, pneumonia and myocardial infarction.

Recent studies describe an association between inflammatory reactions and cardiovascular diseases like arteriosclerosis or latent and chronic persistent infections. As a marker for inflammation, CRP high-sensitive can be used to predict the risk of myocardial infarction and stroke.

The CRP kit provides an easy-to-use assay for monitoring anti-rejection therapies. The ELISA assay was used for years to test hundreds of patients. Its predictive diagnostic value was compared with the gold standard (kidney biopsy).

Indications

- Prognosis factor for myocardial infarction or stroke
- Inflammatory processes

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9710s	PLATE	Micro titre plate with precoated strips	12 x 8 wells
K 9710s	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K 9710s	CONJ	POD antibody, (rabbit-anti-CRP, Peroxidase-labeled)	1 x 150 µl
K 9710s	STD	Calibrators, ready to use* (0; 1.9; 5.6; 16.7; 50; 150 ng/ml)	6 x 1 ml
K 9710s	CTRL	Control, ready-to-use	1 x 1 ml
K 9710s	CTRL	Control, ready-to-use	1 x 1 ml
K 9710s	SAMPLEBUF	Sample buffer, ready to use	2 x 100 ml
K 9710s	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml

Cat. No.	Label	Kit components	Quantity
K 9710s	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

*The CRP calibrators were standardized against WHO standard 470.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C ($\geq 18.2 \text{ M}\Omega \text{ cm}$).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash

buffer). The CONJ is stable at **2–8 °C** until the expiry date stated on the label. **Conjugate (1:101 diluted CONJ) is not stable and cannot be stored.**

- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. PREPARATION OF SAMPLES

We recommend DrySpot-ID (catalogue no DZ9020ID or DZ9021ID) as dried blood spot carrier. The moistened cards are stable for 3 weeks at room temperature.

Place the filter paper containing the dried sample in a 1.5 ml tube. Add 1000 µl SAMPLEBUF and mix well. After incubation for 15 minutes at room temperature (15–30 °C) mix again well. For testing in duplicates, pipette 2 x 100 µl of each prepared sample per well.

7. ASSAY PROCEDURE

Principle of the test

This enzyme immunoassay is a sandwich assay for the determination of CRP in dried blood samples. The wells of the microtitre plate are coated with polyclonal antibodies directed against C-reactive protein. In a first incubation step, after dissolving CRP from the dried blood spot, the CRP in the samples is bound to the coated polyclonal rabbit antibodies (in excess). To remove all unbound substances, a washing step is carried out.

In a second incubation step, a peroxidase-labeled CRP antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stopping solution is then added and the color converts to yellow. The intensity of the yellow color is directly proportional to the concentration of CRP in the sample. A dose response curve of the absorbance (at 450 nm) unit vs. concentration is generated. CRP, present in the patient samples, is determined directly from this calibration curve.

The combination of two specific antibodies in the CRP ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of STD /SAMPLE/CTRL (standards/sample/controls) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper..
2.	Add 100 µl of STD (standard), prepared SAMPLE (sample), and CTRL (controls) into the respective wells.
3.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30 °C) on a horizontal mixer.
4.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper..
5.	Add 100 µl conjugate in each well
6.	Cover the plate tightly and incubate for 1 hour at room temperature (15–30 °C) on a horizontal mixer.
7.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper..
8.	Add 100 µl SUB (TMB substrate) in each well.
9.	Incubate for 10–20 minutes* at room temperature (15–30 °C) in the dark*
10..	Add 50 µl STOP (ELISA stop solution) and mix well
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (595 nm, 630nm, 650 nm or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (595 nm, 630nm, 650 nm or 690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Dried blood samples

The obtained results have to be multiplied with the **dilution factor of 60** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

CRP-Concentration, Serum/plasma (according to Pearson et al., 2003):

< 1 mg/l	low CHD-Risk
1–3 mg/l	medium CHD-Risk
> 3 mg/l	high CHD-Risk

If the CRP concentration is found to be higher than 3 mg/l, a second determination should be made within 2 to 3 weeks. If the CRP-concentration is again high, and other reasons are excluded (acute infection, chronic-inflammatory diseases), the obtained CRP concentration can be used for risk stratification in coronary heart disease (CHD) patients. If the CHD risk is high, the lifestyle should be changed together with medical treatment. These normal ranges should be used as a guideline only.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to CRP. The specificity is calculated in percent, based on the cross-reactivity of these compounds.

• Alpha-1-Antitrypsin	0 %
• Lysozym	0 %
• Albumin	0 %
• Other acute phase proteins	0 %

No cross reactivity with CRP in mouse serum was observed.

Analytical Sensitivity

The sensitivity was set as $B_0 + 2 \text{ SD}$. The zero-standard was measured 18 times.
(n = 18)

Sample	CRP mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.005	0.001	0.921

Spiking Recovery

Two samples with known CRP concentration were spiked with four different amounts of CRP and measured. (n = 4)

Unspiked Sample [ng/ml]	Spike [ng/ml]	CRP expected [ng/ml]	CRP measured [ng/ml]
9.8	37.5	47.3	44.5
9.8	18.8	28.6	27.3
9.8	9.4	19.2	18.2
9.8	4.7	14.5	14.3
9.3	37.5	46.8	48.2
9.3	18.8	28.1	26.3
9.3	9.4	18.7	18.0
9.3	4.7	14.0	13.7

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik CRP ELISA test was calculated from 20 replicate determinations of two samples.

Intra-Assay (n = 20)

Sample	CRP mean value [ng/ml]	Intra-Assay CV [%]
1	23.3	6
2	99.4	5.5

The total precision (inter-assay variation) of the Immundiagnostik CRP ELISA test was calculated from data on 2 samples obtained in 15 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay (n = 14)

Sample	CRP Mean value [ng/ml]	Inter-Assay CV [%]
1	22.1	11.6
2	90.4	13.8

Dilution recovery

Two patient samples were diluted and analyzed. The results are shown below:

Sample	Dilution	Expected [µg/ml]	Measured [µg/ml]	Recovery [%]
A	1:100	2.90	2.88	99.3
	1:200	1.45	1.55	106.8
	1:400	0.73	0.83	113.7
	1:800	0.36	0.39	108.3
	1:1600	0.18	0.18	100.0
B	1:200	10.80	10.80	100.0
	1:400	5.40	5.80	107.4
	1:800	2.70	2.90	107.4
	1:1600	1.35	1.61	119.3
	1:3200	0.68	0.83	122.1
	1:6400	0.33	0.35	106.1

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immunodiagnostik AG along with a written complaint.

15. REFERENCES

1. Koenig, W. et al., 2004. C-reactive protein modulates risk prediction based on the Framingham Score: implications for future risk assessment: results from a large cohort study in southern Germany. *Circulation*, **109**(11), pp.1349–53.

2. Pearson, T.A., 2003. American Heart Association Guide for Improving Cardiovascular Health at the Community Level: A Statement for Public Health Practitioners, Healthcare Providers, and Health Policy Makers From the American Heart Association Expert Panel on Population and Pre. *Circulation*, **107**(4), pp.645–651.

3. Ridker, P.M. et al., 2000. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. *The New England journal of medicine*, **342**(12), pp.836–43.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by