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

Hydroxy-Pyridinium-Crosslinks HPLC Kit

*For the determination of hydroxy-pyridinium crosslinks
in urine*

Gültig ab / Valid from 2015-08-31

KC3201

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE _____ 15

2. INTRODUCTION _____ 15

3. PRINCIPLE OF THE TEST _____ 15

4. MATERIAL SUPPLIED _____ 16

5. MATERIAL REQUIRED BUT NOT SUPPLIED _____ 16

6. STORAGE AND PREPARATION OF REAGENTS _____ 17

7. PRECAUTIONS _____ 17

8. SPECIMEN COLLECTION AND PREPARATION _____ 18

9. ASSAY PROCEDURE _____ 18

Procedural notes _____ 18

Sample hydrolysis _____ 18

Test procedure _____ 18

Chromatographic conditions _____ 19

10. TREATMENT OF THE COLUMN _____ 19

11. RESULTS _____ 19

Typical chromatogram _____ 20

12. LIMITATIONS _____ 20

13. QUALITY CONTROL _____ 20

Expected values _____ 20

14. PERFORMANCE CHARACTERISTICS _____ 21

Precision and reproducibility _____ 21

Linearity _____ 21

Detection limit _____ 21

Recovery _____ 21

15. DISPOSAL _____ 21

16. TROUBLESHOOTING _____ 22

17. REFERENCES _____ 23

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE _____ 23

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of pyridinoline und deoxypyridinoline in urine. This assay is designed for *in vitro* diagnostic use only.

2. INTRODUCTION

The hydroxypyridinium compounds pyridinoline (PYD) and deoxypyridinoline (DPD) are specific constituents of mature skeletal collagens. They are released into the circulation and excreted in the urine. Their measurement in urine is a sensitive index of the extent of ongoing bone resorption. The quantification of collagen crosslinks in urine is achieved by chromatographic techniques.

Clinical applications of hydroxypyridinium markers include numerous metabolic bone disorders such as osteoporosis, primary hyperparathyroidism, Paget's disease of bone, and metastatic bone disease. Urinary pyridinium crosslinks of collagen also show great promise as markers of therapeutic efficacy in bone disorders associated with accelerated bone resorption.

3. PRINCIPLE OF THE TEST

The first step in determining the hydroxy-pyridinium crosslinks is a hydrolysis of the sample followed by a solid phase extraction on cellulose cartridges. 100 µl of the eluate are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 20 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered urine calibrator; the concentration is calculated via integration of the peak areas by the external standard method.

Summary

Besides many other parameters, the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immunoassays with up to six calibrators per test, a one-point calibration is mostly sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher sample numbers can be handled nearly without control.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KC3201LM	MOPHA	Mobile phase	1000 ml
KC3201KA	CAL	Calibrator, lyophilized	3 vials
KC3201HY	HYDREA	Hydrolysis reagent	100 ml
KC3201CS	CELLULOSE	Cellulose	3 x 4 g
KC3201RB	REABUF	Reaction buffer Reagent I Reagent II	3 x 22 ml 3 x 88 ml
KC3201W1	WASH1	Washing solution I Reagent I Reagent II	3 x 280 ml 3 x 560 ml
KC3201W2	WASH2	Washing solution II	50 ml
KC3201EL	ELUSOL	Elution solution	50 ml
KC3201KO	CTRL1 CTRL2	Control 1 and 2; 2.2 ml lyophilized	2 x 3 vials

The HPLC column (KC 3201RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Glass tubes for hydrolysis (dark, heat-sealable)
- Vacuum chamber for solid phase extraction
- HPLC-pump (isocratic) with an Fluorescence-detector
- Superspher RP18-column, 4 µm, 125 × 4 mm
- Cryostat or heating block
- Chromatographic-cartridges with discs or glass wool

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

6. STORAGE AND PREPARATION OF REAGENTS

Preparation of the calibrator and controls

Reconstitute the calibrator (CAL; urine spiked with defined amounts of pyridinoline and desoxypyridinoline) in 5.5 ml ultra pure water. Divide the calibrator (CAL) in several aliquots and store them at -20 °C. Avoid repeated freeze-thaw cycles. The contents of pyridinium and desoxypyridinium crosslinks slightly change from lot to lot. The exact content is given on the label.

Reconstitute the controls (CTRL1, CTRL2) in 2.2 ml ultra pure water.

All other test reagents are stable at 20–25 °C, up to the date of expiry stated on the label.

Preparation of the reaction buffer

Add the total content of reaction buffer (REABUF) reagent I to the bottle of reaction buffer (REABUF) reagent II and incubate the mixture for 12 hours at room temperature. The freshly prepared reaction buffer is stable for 3 months.

Preparation of washing solution I

Add the total content of washing solution I (WASH1) reagent I to the bottle of washing solution (WASH1) reagent II and incubate the mixture for 12 hours at room temperature. The freshly prepared washing solution I is stable for 3 months.

Preparation of the cellulose slurry

Add 40 ml of washing solution I to 4 g cellulose (CELLULOSE) and mix well. The freshly prepared cellulose slurry is stable for 1 month.

7. PRECAUTIONS

- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should nevertheless be handled as if potentially infectious.
- The hydrolysis reagent (HYDREA), reaction buffer (REABUF) and washing solution I (WASH1) contain acid. Although diluted, they still must be handled with care. They can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapors and avoid inhalation.

- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Urine is suited for this test system. A two hour urine sampling between 7⁰⁰ and 10⁰⁰am correlates well with a 24 h collection. Samples are stable for 24 h at room temperature and up to one week at 2–8 °C. Longer storage should be at -20 °C.

9. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- The assay should always be performed according the enclosed manual.

Sample hydrolysis

Patient urine, calibrator (CAL) or controls (CTRL1, CTRL2) should be diluted 1+1 (1 ml + 1 ml) with hydrolysis reagent (HYDREA) and hydrolysed for 4 hours at 120 °C.

Caution: HYDREA consists of concentrated hydrochloride acid.

Test procedure

1.	Fill up the column with 1 ml cellulose slurry. Mix the slurry carefully so that the column is packed homogeneously.
2.	Soak 5 ml washing solution I through the column by vacuum.
3.	Mix 0.5 ml washing solution I, 0.5 ml hydrolyzed cold urine and 3 ml reaction buffer. Shake very well, pipet onto the column and let it being soaked through by vacuum.
4.	Rinse the column with 2 x 8 ml washing solution I by vacuum.
5.	Dry the column by vacuum for 10 min.

6.	Rinse the column with 0.5 ml washing solution II (WASH2) by vacuum.
7.	Dry the column by vacuum for 10 min.
8.	Elute the sample with 0.5 ml elution solution (ELUSOL) by vacuum.
9.	Centrifuge the sample and inject 100 µl into the HPLC system.

Chromatographic conditions

Column material: Merck Supersher, 4 µm

Column dimension: 125 × 4 mm

Flow rate: 1–1.5 ml/min

Temperature: 30 °C

Fluorescence de- Excitation 290 nm
tection: Emission 400 nm

Injection volume: 100 µl

Running time / chromatogram: 20 minutes

10. TREATMENT OF THE COLUMN


After analysis the column should be flushed with 30 ml ultra pure water (1 ml/min) and stored in 50 % methanol in ultra pure water (approx. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 30 ml eluent.

11. RESULTS

$$\frac{\text{Peak height sample} \times \text{calibrator concentration (pmol/ml)}}{\text{Peak height calibrator}} = \text{sample conc. (pmol/ml)}$$

Tip: Alternatively, the peak area instead of the peak height can be used for quantification.

Typical chromatogram

12. LIMITATIONS

Don't use serum, plasma or whole blood.

13. QUALITY CONTROL

Expected values

38.4 ± 18.4 nmol PYD/mmol creatinine

10.3 ± 5.3 nmol DPD/mmol creatinine

We recommend each laboratory to establish its own reference range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or urine pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV

PYD	4.5 % (958 pmol/ml)	[n = 10]
DPD	6.2 % (238 pmol/ml)	[n = 10]

Inter-Assay CV

PYD	8.9 % (1051 pmol/ml)	[n = 10]
PYD	7.9 % (177 pmol/ml)	[n = 10]
DPD	8.6 % (254 pmol/ml)	[n = 10]
DPD	8.2 % (45 pmol/ml)	[n = 10]

Linearity

PYD	up to 2500 pmol/ml
DPD	up to 1000 pmol/ml

Detection limit

PYD	7.3 pmol/ml
PYD	7.3 pmol/ml

Recovery

PYD	97.4 %
DPD	94.4 %

15. DISPOSAL

The mobile phase (MOPHA), washing solution I and II (WASH1, WASH2), reaction buffer (REABUF) and elution solution (ELUSOL) must be disposed as non-halogenated solvent. The hydrolysis solution (HYDREA) can be neutralized with NaOH. If the pH value is neutral, it can be disposed as salt solution.

Important: Reaction will produce extreme heat, be careful!

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Doublepeaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not yet reach working temperature	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system

Problem	Possible reason	Solution
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

17. REFERENCES

1. Seibel, M.J., Robins, S.P. & Bilezikian, J.P., 1992. Urinary pyridinium crosslinks of collagen: Specific markers of bone resorption in metabolic bone disease. *Trends in Endocrinology and Metabolism*, **3**(7), pp.263–270.
2. Robins, S.P. et al., 1991. Evaluation of urinary hydroxypyridinium crosslink measurements as resorption markers in metabolic bone diseases. *European journal of clinical investigation*, **21**(3), pp.310–315.
3. Seibel, M.J. et al., 1992. Urinary hydroxypyridinium cross-links of collagen in primary hyperparathyroidism. *The Journal of clinical endocrinology and metabolism*, **74**(3), pp.481–6.
4. Zittermann, A. et al., 2002. Oral contraceptives moderately effect bone resorption markers and serum-soluble interleukin-6 receptor concentrations. *Calcified tissue international*, **70**(1), pp.11–21.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the

test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by