

IDK[®] Homocystein DBS LC-MS/MS

***Zur Bestimmung von Homocystein in Trockenblutspots
For the determination of homocysteine in dried blood spots***

Gültig ab 2017-02-17

KM7000DBS

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Inhalt

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IDK[®] Homocystein DBS



LC-MS/MS

For the determination of homocysteine in dried blood spots

Valid from 2017-02-17

KM7000DBS

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	17
7. STORAGE AND PREPARATION OF SAMPLES	18
8. TEST PROCEDURE	18
9. CHROMATOGRAPHIC CONDITIONS	19
<i>HPLC method</i>	19
10. MS/MS METHOD	20
<i>MRM transitions (m/z) for Waters Micromass Quattro Ultima</i>	20
<i>MRM transitions (m/z) for Waters XEVO TQ MS</i>	20
11. CALCULATION	20
12. CHROMATOGRAM EXAMPLES	21
13. QUALITY CONTROL	21
<i>Reference range (whole blood)</i>	22
14. PERFORMANCE CHARACTERISTICS	22
<i>Precision and reproducibility</i>	22
<i>Detection limit</i>	22
15. PRECAUTIONS	22
16. TECHNICAL HINTS	23
17. DISPOSAL	23
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
19. REFERENCES	24

1. INTENDED USE

This *dried blood spot* (DBS) LC-MS/MS application is intended for the quantitative determination of homocysteine in dried blood spots. For *in vitro* diagnostic use only.

2. INTRODUCTION

Homocysteine occurs in blood as well as free molecule as bound to albumin. Following a distinct metabolic pathway, methionin is changed to homocystein, which then is further metabolized by cystathionine beta synthase (cofactor vitamin B6) resulting in cystein which is then degraded. A deficiency of the enzyme results in rising concentrations of methionin and homocysteine and a decreasing concentration of cystine. This leads to homocystinuria which causes psychomotor disabilities if untreated.

It has been shown that homocysteine is an important risk factor for coronary heart disease. High values can be found in samples of patients with a genetic defect which results in a failure in homocysteine metabolism. It has been published that this is the cause of 40% of all heart attacks in the USA. Furthermore, the risk of stroke increases proportionally to the concentration of homocystein in blood.

Vitamin B6, B12 and folic acid can be used for reducing high plasma levels of homocystein effectively and without known side-effects.

3. PRINCIPLE OF THE TEST

In the first step, the homocystein in the dried blood spots is reduced and extracted. Therefore reduction solution including an internal standard is added. In the second step, the sample is purified by filter centrifugation plates, diluted and injected.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KM7000LA	MOPHA A	Mobile phase A	1000 ml
KM7000LB	MOPHA B	Mobile phase B	1000 ml
KM7000ST	STD	Standard (lyoph. 200 µl); see specification for concentration	4 vials
KM7000KO	CTRL1 CTRL2	Dried blood spots on dried blood spot carrier; see specification for concentrations	2 x 8
KM7000RL	REDSOL	Reduction solution (lyoph. 600 µl)	4 vials
KM7000VL	DIL	Dilution reagent	22 ml
KM7000AC	ACTSOL	Activation solution	2.5 ml
KM7000IS	INT STD	Internal standard	4.5 ml
KM7000SO	SOL	Standard dilution solution	5 ml

The HPLC separation column (KM7000RP), the elution plates (KM7000PL) and the 96-well plates (KM7000MSP) can be ordered separately from Immundiagnostik AG. The complete kit well as all individual components can be ordered separately. Please ask for our single component price list.

For the preparation of the *IDK*® homocysteine LC-MS/MS application, tuning is necessary for finding the optimal LC-MS/MS settings. The pure substances for the tuning (KM7000TU [5 ml], KM7000TS [2 ml]) can also be ordered at Immundiagnostik AG.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipettors and 10–1000 µl tips
- Multi-channel pipet
- Centrifuge, 3000 g, for elution plates and 96-well plates
- Vortex
- HPLC column, Atlantis™ HILIC Silica 2.1 x 50 mm, 3 µm (KM7000RP)
- Pre-column, Atlantis™ HILIC Silica 3 µm, 2.1 x 10 mm Guard Column (KM-7000VS)
- Guard column holder (KM7000VH)
- LC-MS/MS machine
- Elution plate KM7000PL, multiscreen filter plates (Merck #MSGVN2210)
- 96-well plate KM7000MSP (Greiner #655201)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.**
- Reconstitute the **reduction solution** (REDSOL) with **600 µl ultra pure water**
- **The mobile phases** (MOPHA A, MOPHA B) and **dilution reagent** (DIL) must be prepared before use by adding **0.1 % activation solution** (ACTSOL), e. g.

500 ml MOPHA + 500 µl ACTSOL

22 ml DIL + 22 µl ACTSOL

Reagents containing ACTSOL can be stored for 4 weeks at 2–8 °C.

ATTENTION: The **activation solution** must be added under a **fume hood**. All vessels used must be absolutely clean, free of detergents and preferably made of LC-MS/MS-suitable glass.

- **Controls** (CTRL) are delivered in the kit as dried blood spots on filter paper. The concentration is given on the respective specification. Controls can be stored at **2–8 °C** until the expiry date stated on the label.
- **Preparation of the standard:** The supplied standard has to be diluted to a **standard dilution series**. Reconstitution details and a recommendation for the dilution are given in the **specification data sheet**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

7. STORAGE AND PREPARATION OF SAMPLES

- Suitable samples are dried capillary blood spots on DrySpot-ID filter paper. (catalogue no. DZ9020ID) DrySpot-ID is the only carrier for dried blood samples that can be used. The dried blood spots are stable for 3 months at 2–8 °C.
- For each sample, one dried blood spot is punched out of the card (two for duplicates) and put into the corresponding well of the plate. The size of the dried blood spots is 1/8" / 3.2 mm and equates to 3 µl whole blood.
- Not needed elution plate wells can be stored covered.

8. TEST PROCEDURE

We recommend to carry out the tests in duplicate.

1.	Put elution plate on 96-well plate.
2.	Put dried blood spots from samples and controls (CTRL) in elution plate wells , as well as 30 µl freshly diluted standards (STD).
3.	Mix internal standard (INT STD) and reduction solution (REDSOL) at the ratio of 2:1 , e.g. 1000 µl INT STD + 500 µl REDSOL.
4.	Add 60 µl of the mixture to each well.
5.	Shake for 30 min at room temperature in the dark.
6.	Centrifuge for 2 min at 3000 g.
7.	Remove elution plate and discard it if used up.
8.	Add 200 µl dilution reagent (DIL) to the eluate.
9.	Cover 96-well plate with aluminium foil and press it on tightly to prevent evaporation.
10.	Put 96-well plate in LC-MS/MS-System.
11.	Inject 3 µl.

9. CHROMATOGRAPHIC CONDITIONS

HPLC method

Column: Atlantis™ HILIC Silica 2.1 x 50 mm, 3 µm (Waters 186002011)

Flow: 300–350 µl/min

Column temperature: 30°C

Injection volume: 3 µl

Gradient:

Time (min)	% A	% B	Flow rate [µl/min]
0.0	95	5	300
0.5	95	5	300
4.5	70	30	350
5.0	60	40	350
7.0	95	5	350

We recommend to use a guard column/filter to extend column's life. The guard column cartridge (KM7000VS) and guard column holder (KM7000VH) can be ordered separately from Immundiagnostik.

In addition, we recommend to use a shift valve to preserve the mass spectrometer; the first 2.5 min of each run can then go to the waste.

After the analysis, the column should be washed with acetonitril at 300 µl/min for about 30 min.

10. MS/MS METHOD

The parameters for the MS/MS method are listed as examples for two machines in the following table.

Column A shows the parameters for a **Waters Micromass Quattro Ultima Tandem mass spectrometer**.

Column B shows the parameters for a **Waters XEVO TQ MS**.

Parameter	A	B
Mode:	MRM	MRM
Polarity:	ESI+	ESI+
Capillary (kV):	3,5	3,0
Cone (V):	20	20
Extracor (V):	3	3
RF lens (V):	3	–
Source temperature (°C):	100	150
Desolvation temperature (°C):	250	600
Cone gas flow (L/Hr):	120	20
Desolvation gas flow (L/Hr):	700	100

MRM transitions (m/z) for Waters Micromass Quattro Ultima

Homocysteine	136 > 90	Collision energy: 11
DL-d ₄ -homocysteine	140 > 94	Collision energy: 11

MRM transitions (m/z) for Waters XEVO TQ MS


Homocysteine	136 > 90	Collision energy: 12
DL-d ₄ -homocysteine	140 > 94	Collision energy: 10

11. CALCULATION


Linear regression in consideration of the internal standard can be used as model for evaluation of the results. The standard concentration points are connected by a straight line. The samples can be calculated using the obtained line.

12. CHROMATOGRAM EXAMPLES

Sample

Calibrator

13. QUALITY CONTROL

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range (whole blood)

< 8 µmol/l	no increased homocystein-dependent risk
8–10 µmol/l	borderline values. Negative effects, for example by influence of coagulation, already possible. Increased risk not to be ruled out, therapy recommended.
10–15 µmol/l	increased homocysteine level with aggressive vascular effects; clearly increased risk, e.g. for vasoconstriction in the carotis, shown by studies. Therapy necessary.
> 15 µmol/l	clearly increased homocysteine level with clearly increased risk for atherosclerosis in the heart, brain, peripheral arteries, increased risk for thromboembolism; therapy absolutely necessary. Duplication of risk for Alzheimer's and Parkinson's disease.

Based on Immundiagnostik studies of samples of apparently healthy persons (n = 9), a normal range of 2.8–7.9 µmol/l was estimated.

We recommend each laboratory to establish its own reference range.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV

4.8 % (19.8 µmol/l) [n = 10]

5.1 % (25.5 µmol/l) [n = 10]

Inter-Assay CV

9.6 % (8.99 µmol/l) [n = 5]

Detection limit

Detection limit of homocysteine: 0.42 µmol/l

Please note that the determination of the detection limit depends not only on the application method but also on the instrument.

15. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- The activation solution consists of acid. It must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

16. TECHNICAL HINTS

- Quality control samples should be analyzed with each run.
- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date stated on kit label.
- The assay should always be performed according the enclosed manual.

17. DISPOSAL

- The mobile phases (MOPHA A, MOPHA B), standard dilution solution (SOL), dilution reagent (DIL) and activation solution (ACTSOL) must be disposed as non-halogenated solvents.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK*® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

19. REFERENCES

<http://www.deutsche-gefaessliga.de/homocystein.html>

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by

Attention